


K5120 KBBC Medium

Crop: **Bean, Leek, Pea, Tomato**
 Disease: **Bacterial brown spot (bean)**
 Pathogen: ***Pseudomonas syringae* pv. *syringae***
Pseudomonas syringae* pv. *porri
Pseudomonas syringae* pv. *pisii
Pseudomonas syringae* pv. *tomato



Pseudomonas syringae pv. *syringae* (*Pss*) is the causal organism of bacterial brown spot of beans. This bacterium is seed borne and therefore its detection on seeds is important. KBBC medium is a rather selective medium to detect *Pss* on seeds of beans. This medium is based on King's B Medium (K5165), however in KBBC Medium boric acid (1.5 g/liter), cephalixin and nystatin are added. Nystatin is used to control fungi. As an alternative, cycloheximide, a more potent fungicide, can be used. KBBC is much more selective than MSP (M5167) and in general the recovery of *Pss* is smaller on KBBC than on MSP. *Pspha*, unlike *Pss*, will not grow on KBBC. Therefore, the chance of detection of *Pss* is higher when both complementary media are used. Detection of *Pss* is performed by the dilution plating of bacterial extract on KBBC and MSP. Then *Pss*-suspected isolates are transferred to KB medium. Finally, the identification of suspected colonies can be performed by a pathogenicity assay or PCR. Colonies of *Pss* on KBBC are 3-4 mm in diameter, flat, circular, translucent, creamy white and show blue fluorescence under UV light. This medium can also be used for the detection of seed borne *Ps porri*, *Ps pisi* and *Ps tomato* on seed of resp. leek, pea and tomato.

COMPOSITION OF MEDIA K5120: KBBC MEDIUM

COMPOUND	GRAM/LITER
Agar	15.0
Di-potassium hydrogen phosphate (K ₂ HPO ₄)	1.5
Boric acid (H ₃ BO ₃)	1.5
Magnesium sulphate anhydrous (MgSO ₄ anhydrous)	0.73
Proteose Peptone	20.0

METHOD

- Dissolve 38.7 grams of ingredients in distilled water and adjust volume to 970 ml.
- Add 30 ml glycerol (50%) and mix.
- Adjust pH to 7.2.
- Autoclave the solutions (121 °C, 15 psi, 15 minutes).
- Prepare sterile antibiotic solutions and add the following amounts per liter medium:
80 mg cephalixin monohydrate (C0110)
35 mg nystatin (N0138) or 100 mg cycloheximide (C0176)
- Allow medium to cool down to ca. 45 °C – 50 °C and add antibiotics to the solution.
- Mix gently to avoid air bubbles and pour plates (15-20 ml per 9.0 cm plate).

Reference:

Mohan, S.K. and Schaad, N.W. 1987. An improved agar plating assay for detecting *Pseudomonas syringae* pv. *syringae* and *Pseudomonas syringae* pv. *phaseolicola* in contaminated bean seed. *Phytopathology* 77: 1390-1395.

K5120 KBBC MEDIUM

K5120.1000 1 kg

For prepared and ready to use plates of this medium contact:
 Tritium Microbiologie Tel : 040-2051615
 Rooijakkersstraat 6 Fax : 040-2051395
 5652 BB Eindhoven Email : info@tritium-microbiologie.nl
 The Netherlands