

FICHE TECHNIQUE

A0101 : AMOXYCILLIN TRIHYDRATE

CAS number 61336-70-7

Molecular weight $C_{16}H_{19}N_3O_5S \cdot 3H_2O = 419.5$

Assay > 95 %

Storage condition

- Store at room temperature
- Soluble in water

Hazard statements

- **H317** May cause an allergic skin reaction
- **H334** May cause allergy or asthma symptoms or breathing difficulties if inhaled
- **H402** Harmful to aquatic life

Precautionary statements

- **P261** Avoid breathing dust/fume/gas/mist/vapours/spray
- **P273** Avoid release to the environment
- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P302** IF ON SKIN:
- **P311** Call a POISON CENTER or doctor/physician
- **P342** If experiencing respiratory symptoms:
- **P352** Wash with soap and water

A0101.0010/ 10 g

A0101.0025/ 25 g