

FICHE TECHNIQUE

A0164: APRAMYCIN SULPHATE

CAS number 65710-07-8

Molecular weight $C_{21}H_{41}N_5O_{11} \cdot 2.5H_2SO_4 = 784.8$

Assay > 50 % (base)

Storage condition

- Store at 2-8 °C
- Soluble in water, protect from moisture

Hazard statements

- **H315** Causes skin irritation
- **H319** Causes serious eye irritation
- **H360** May damage fertility or the unborn child

Precautionary statements

- **P201** Obtain special instructions before use
- **P305** IF IN EYES:
- **P308** IF exposed or concerned:
- **P313** Get medical advice/attention
- **P351** Rinse continuously with water for several minutes

A0164.0005/ 5 g

A0164.0010/ 10 g

A0164.0025/ 25 g