

FICHE TECHNIQUE



A0201: ANDERSON'S RHODODENDRON

Micro Elements

mg/l

μM

CoCl₂·6H₂O

0.025

0.11

CuSO₄·5H₂O

0.025

0.10

FeNaEDTA

73.40

200.00

H₃BO₃

6.20

100.27

KI

0.30

1.81

MnSO₄·H₂O

16.90

100.00

Na₂MoO₄·2H₂O

0.25

1.03

ZnSO₄·7H₂O

8.60

29.91

Macro Elements

mg/l

mM

CaCl₂

332.02

2.99

KNO₃

480.00

4.75

MgSO₄

180.54

1.50

NaH₂PO₄

330.60

2.75

NH₄NO₃

400.00

5.00

Anderson W.C., Tissue culture propagation of Red and Black Raspberries, *Rubus Idaeus* and *R. Occidentalis* Act. Hort., 112, 13 (1980).

Hazard statements

- **H315** Causes skin irritation
- **H319** Causes serious eye irritation
- **H335** May cause respiratory irritation

Precautionary statements

- **P261** Avoid breathing dust/fume/gas/mist/vapours/spray
- **P302** IF ON SKIN:
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes
- **P352** Wash with soap and water

A0201.0001/ 1 l

A0201.0010/ 10 l

A0201.0050/ 50 l

