Safety Data Sheet

according to Regulation (EU) 2015/830 Date of issue: 20/06/2012 Revision date 12/07/2018 A0202

Version: 2.0

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier	
Product form	: Mixture
Trade name	: Anderson's Rhododendron
	(Micro and Macro elements incl. vitamins)
Product code	: A0202
Product group	: Raw material

1.2. Relevant identified uses of the substance or mixture and uses advised against

1.2.1. Relevant identified uses

Main use category	: Professional use
Industrial/Professional use spec	: For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes.

1.2.2. Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

Manufacturer

Duchefa Biochemie B.V. A. Hofmanweg 71 2031 BH Haarlem - The Netherlands T +31(0)23-5319093 - F +31(0)23-5318027 info@duchefa.nl

1.4. Emergency telephone number

Emergency number

: Supplier contact information: +31(0)23-5319093 (M-F 09:00-17:00)

+31(0)6-30109355 (outside office hours)

Organisation/Company	Address	Comment
World Health Organization world directory of poison centres	http://apps.who.int/poisoncentres/	Consult website for a local poison centre

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Skin corrosion/irritation, Category 2	H315
Serious eye damage/eye irritation, Category 2	H319
Specific target organ toxicity — Single exposure, Category 3,	H335
Respiratory tract irritation	

Full text of H statements : see section 16

Adverse physicochemical, human health and environmental effects

No additional information available

Safety Data Sheet

according to Regulation (EU) 2015/830

2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

Hazard pictograms (CLP)	GHS07
Signal word (CLP)	: Warning
Hazard statements (CLP)	 H315 - Causes skin irritation. H319 - Causes serious eye irritation. H335 - May cause respiratory irritation.
Precautionary statements (CLP)	 P261 - Avoid breathing dust. P302+P352 - IF ON SKIN: Wash with plenty of water. P305+P351 - IF IN EYES: Rinse cautiously with water for several minutes.
Extra phrases	: Based on research by TNO in Rijswijk (The Netherlands), commissioned by Duchefa Biochemie B.V. in Haarlem, the medium has no oxidising or explosive properties. As such the substance is not classified as oxidizing (H272, GHS03).

2.3. Other hazards

No additional information available

SECTION 3: Composition/information on ingredients

3.1. Substances

Not applicable

3.2. Mixtures

Name	Product identifier	%	Classification according to Regulation (EC) No. 1272/2008 [CLP]
Potassium nitrate	(CAS-No.) 7757-79-1 (EC-No.) 231-818-8	23,889	Ox. Sol. 2, H272
Ammonium nitrate	(CAS-No.) 6484-52-2 (EC-No.) 229-347-8	19,908	Ox. Sol. 3, H272
Calcium chloride	(CAS-No.) 10043-52-4 (EC-No.) 233-140-8 (EC Index-No.) 017-013-00-2	16,524	Eye Irrit. 2, H319
Sodium dihydrogenorthophosphate	(CAS-No.) 7558-80-7 (EC-No.) 231-449-2	16,454	Not classified
Magnesium sulfate anhydrous	(CAS-No.) 7487-88-9 (EC-No.) 231-298-2	8,985	Not classified
Myo-Inositol	(CAS-No.) 87-89-8 (EC-No.) 201-781-2	4,977	Not classified
Adenine hemisulphate	(CAS-No.) 321-30-2 (EC-No.) 206-286-5	3,982	Acute Tox. 4 (Oral), H302
Ethylenediaminetetraacetate ferric sodium	(CAS-No.) 15708-41-5 (EC-No.) 239-802-2	3,653	Not classified

Safety Data Sheet

according to Regulation (EU) 2015/830

Manganese sulfate monohydrate	(CAS-No.) 10034-96-5 (EC-No.) 232-089-9 (EC Index-No.) 025-003-00-4	0,841	Eye Dam. 1, H318 STOT RE 2, H373 Aquatic Chronic 2, H411
Zinc sulphate	(CAS-No.) 7446-20-0 (EC-No.) 231-793-3 (EC Index-No.) 030-006-00-9	0,428	Acute Tox. 4 (Oral), H302 Eye Dam. 1, H318 Aquatic Acute 1, H400 Aquatic Chronic 1, H410
Boric acid substance listed as REACH Candidate	(CAS-No.) 10043-35-3 (EC-No.) 233-139-2 (EC Index-No.) 005-007-00-2	0,309	Repr. 1B, H360FD
Thiamine hydrochloride	(CAS-No.) 67-03-8 (EC-No.) 200-641-8	0,02	Not classified
Potassium iodide	(CAS-No.) 7681-11-0 (EC-No.) 231-659-4	0,0149	Acute Tox. 4 (Oral), H302 Skin Irrit. 2, H315 Eye Irrit. 2, H319
Sodium molybdate dihydrate	(CAS-No.) 10102-40-6 (EC-No.) 231-551-7	0,0124	Not classified
Cobalt(II) chloride substance listed as REACH Candidate (Cobalt dichloride)	(CAS-No.) 7646-79-9 (EC-No.) 231-589-4 (EC Index-No.) 027-004-00-5	0,0012	Acute Tox. 4 (Oral), H302 Resp. Sens. 1, H334 Skin Sens. 1, H317 Muta. 2, H341 Carc. 1B, H350i Repr. 1B, H360F Aquatic Acute 1, H400 Aquatic Chronic 1, H410
Copper sulphate	(CAS-No.) 7758-99-8 (EC-No.) 231-847-6 (EC Index-No.) 029-004-00-0	0,0012	Acute Tox. 4 (Oral), H302 Skin Irrit. 2, H315 Eye Irrit. 2, H319 Aquatic Acute 1, H400 Aquatic Chronic 1, H410

Specific concentration limits:

Name	Product identifier	Specific concentration limits
Boric acid	(CAS-No.) 10043-35-3 (EC-No.) 233-139-2 (EC Index-No.) 005-007-00-2	(C >= 5,5) Repr. 1B, H360FD
Cobalt(II) chloride	(CAS-No.) 7646-79-9 (EC-No.) 231-589-4 (EC Index-No.) 027-004-00-5	(C >= 0,01) Carc. 1B, H350i

Full text of H-statements: see section 16

SECTION 4: First aid measures

4.1. Description of first aid measures

: Seek medical attention if ill effect develops.
: Remove victim to fresh air.
: Wash skin with mild soap and water.
: Rinse with water.
: Rinse mouth.

Safety Data Sheet

4.2. Most important symptoms and effects, both acute and delayed

Symptoms/effects after skin contact	: Causes skin irritation.
Symptoms/effects after eye contact	: Redness, pain.

4.3. Indication of any immediate medical attention and special treatment needed

No additional information available

SECT	FION 5: Firefighting mea	asures
5.1. Suitable	Extinguishing media e extinguishing media	 Alcohol resistant foam. Dry chemical powder. Carbon dioxide (CO2). Water spray.
5.2. Hazard case of	fire	 the substance or mixture Under fire conditions, hazardous fumes will be present: COx NOx SOx.
-	Advice for firefighters ating instructions ion during firefighting	Prevent fire fighting water from entering the environment.Wear proper protective equipment.
SECT	FION 6: Accidental relea	ise measures
	Personal precautions, protec al measures For non-emergency personne	tive equipment and emergency procedures : Avoid raising powdered materials into airborne dust.
Emerge	ency procedures	: Wear suitable protective clothing.
	For emergency responders itional information available	
6.2. Preven	Environmental precautions t entry to sewers and public wate	ers.
6.3. Method	Methods and material for con Is for cleaning up	tainment and cleaning up : Sweep up dry powder and dispose properly.
6.4. No add	Reference to other sections itional information available	
SECT	TION 7: Handling and st	orage
7.1.	Precautions for safe handling	3

Precautions for safe handling : Avoid dust formation. Handle in accordance with good industrial hygiene and safety procedures.

Safety Data Sheet

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions

: Store at room temperature Store in dry, well-ventilated area Hygroscopic Keep container tightly closed and dry.

7.3. Specific end use(s)

For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

No additional information available

8.2. Exposure controls

Hand protection:

Туре	Material	Permeation	Thickness (mm)	Standard
Gloves	Nitrile rubber (NBR)	6 (> 480 minutes)	0,11	EN 374

Eye protection:

Safety glasses (to European standard EN 166 or equivalent)

Skin and body protection:

In case of possible repeated skin contact wear protective clothing

Respiratory protection:

Where excessive dust may result, wear approved mask. Type P2 (EN 143)

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state Appearance Colour	:Solid :Powder. :White to slightly yellow.
Odour Odour threshold	: Characteristic. Weak. : No data available
pH	: No data available
Relative evaporation rate (butylacetate=1)	: No data available
Melting point	: No data available
Freezing point	: No data available
Boiling point	: No data available
Flash point	: No data available
Auto-ignition temperature	: No data available
Decomposition temperature	: No data available
Flammability (solid, gas)	: No data available
Vapour pressure	: No data available
Relative vapour density at 20 °C	: No data available

Safety Data Sheet

according to Regulation (EU) 2015/830

Relative density Solubility	: No data available : Readily soluble in water.
Log Pow Viscosity, kinematic Viscosity, dynamic Explosive properties Oxidising properties	 No data available
Explosive limits	: No data available

9.2. Other information

No additional information available

SECTION 10: Stability and reactivity

10.1. Reactivity

Stable under normal conditions of storage, handling and use.

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

No additional information available

10.4. Conditions to avoid

Moisture.

10.5. Incompatible materials

Strong oxidizers.

10.6. Hazardous decomposition products

Thermal decomposition generates :

- COx
- NOx
- SOx

- Phosphorus oxides.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity	: Not classified
Skin corrosion/irritation Serious eye damage/irritation Respiratory or skin sensitisation Germ cell mutagenicity Carcinogenicity	 Causes skin irritation. Causes serious eye irritation. Not classified Not classified Not classified Not classified
Reproductive toxicity STOT-single exposure STOT-repeated exposure	Not classifiedMay cause respiratory irritation.Not classified

Safety Data Sheet

according to Regulation (EU) 2015/830

Aspiration hazard

: Not classified

SECTION 12: Ecological information

12.1. Toxicity

No additional information available

12.2. Persistence and degradability

No additional information available

12.3. Bioaccumulative potential

No additional information available

12.4. Mobility in soil

No additional information available

12.5. Results of PBT and vPvB assessment

Anderson's Rhododendron (Micro and Macro elements incl. vitamins)	
Component	
Boric acid (10043-35-3)	This substance/mixture does not meet the PBT criteria of REACH regulation, annex XIII This substance/mixture does not meet the vPvB criteria of REACH regulation, annex XIII

12.6. Other adverse effects

Additional information

Waste treatment methods

: Prevent entry to sewers and public waters. Avoid release to the environment

SECTION 13: Disposal considerations

13.1. Waste treatment methods

: Dispose in a safe manner in accordance with local/national regulations. Avoid release to the environment.

SECTION 14: Transport information

In accordance with ADR / IATA / IMDG

ADR	IMDG	IATA	
14.1. UN number	14.1. UN number		
Not regulated	Not regulated	Not regulated	
14.2. UN proper shipping name			
Not regulated	Not regulated	Not regulated	
Not regulated	Not regulated	Not regulated	
14.3. Transport hazard class(es)			
Not regulated	Not regulated	Not regulated	
Not regulated	Not regulated	Not regulated	
14.4. Packing group			
Not regulated	Not regulated	Not regulated	
14.5. Environmental hazards			
Not regulated	Not regulated	Not regulated	

Safety Data Sheet

according to Regulation (EU) 2015/830

ADR	IMDG	ΙΑΤΑ
No supplementary information available		

14.6. Special precautions for user

Overland transport

Not regulated

- Transport by sea Not regulated

- Air transport

Not regulated

14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

15.1.1. EU-Regulations

Contains no REACH substances with Annex XVII restrictions Contains a substance on the REACH candidate list in concentration ≥ 0.1% or with a lower specific limit: Boric acid (EC 233-139-2, CAS 10043-35-3)

Contains no REACH Annex XIV substances

15.1.2. National regulations

Ensure all national/local regulations are observed.

Germany

Reference to AwSV 12th Ordinance Implementing the Federal Immission Control Act -	 Water hazard class (WGK) 3, severe hazard to water (Classification according to AwSV, Annex 1) Is not subject of the 12. BlmSchV
12.BImSchV	(Hazardous Incident Ordinance)
Netherlands	
SZW-lijst van kankerverwekkende stoffen	: Manganese sulfate monohydrate,Cobalt(II) chloride are listed
SZW-lijst van mutagene stoffen	: Manganese sulfate monohydrate is listed
NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Borstvoeding	: None of the components are listed
NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Vruchtbaarheid	: Boric acid,Cobalt(II) chloride are listed
NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Ontwikkeling	: Boric acid is listed
Denmark	

Safety Data Sheet

according to Regulation (EU) 2015/830

Recommendations Danish Regulation

: Young people below the age of 18 years are not allowed to use the product

Pregnant/breastfeeding women working with the product must not be in direct contact with the product

The requirements from the Danish Working Environment Authorities regarding work with carcinogens must be followed during use and disposal

15.2. Chemical safety assessment

No additional information available

SECTION 16: Other information

Indication of changes:

2.2	Hazard statements (CLP)	Added	H315;H319;H335
2.2	Signal word (CLP)	Added	Warning
2.2	Hazard pictograms (CLP)	Added	GHS07
2.2	Precautionary statements (CLP)	Added	P261;P302+P352;P305+P351

Abbreviations and acronyms:

	, ,
ATE	Acute Toxicity Estimate
ADR	European Agreement concerning the International Carriage of Dangerous Goods by Road
BCF	Bioconcentration factor
CLP	Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008
DPD	Dangerous Preparations Directive 1999/45/EC
DSD	Dangerous Substances Directive 67/548/EEC
IATA	International Air Transport Association
IMDG	International Maritime Dangerous Goods
LC50	Median lethal concentration
LD50	Median lethal dose
LOAEL	Lowest Observed Adverse Effect Level
NOAEC	No-Observed Adverse Effect Concentration
PBT	Persistent Bioaccumulative Toxic
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006
SDS	Safety Data Sheet

Data sources

: REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006. TNO (Netherlands Organisation for Applied Scientific Research).

Full text of H- and EUH-statements:

Acute Tox. 4 (Oral)	Acute toxicity (oral), Category 4
Aquatic Acute 1	Hazardous to the aquatic environment — Acute Hazard, Category 1
Aquatic Chronic 1	Hazardous to the aquatic environment — Chronic Hazard, Category 1
Aquatic Chronic 2	Hazardous to the aquatic environment — Chronic Hazard, Category 2
Carc. 1B	Carcinogenicity (inhalation) Category 1B
Eye Dam. 1	Serious eye damage/eye irritation, Category 1
5	

12/07/2018 Duchefa Biochemie B.V.

Safety Data Sheet

according to Regulation (EU) 2015/830

Eye Irrit. 2	Serious eye damage/eye irritation, Category 2
Muta. 2	Germ cell mutagenicity, Category 2
Ox. Sol. 2	Oxidising Solids, Category 2
Ox. Sol. 3	Oxidising Solids, Category 3
Repr. 1B	Reproductive toxicity, Category 1B
Repr. 1B	Reproductive toxicity, Category 1B
Resp. Sens. 1	Respiratory sensitisation, Category 1
Skin Irrit. 2	Skin corrosion/irritation, Category 2
Skin Sens. 1	Skin sensitisation, Category 1
STOT RE 2	Specific target organ toxicity — Repeated exposure, Category 2
STOT SE 3	Specific target organ toxicity — Single exposure, Category 3, Respiratory tract irritation
H272	May intensify fire; oxidiser.
H302	Harmful if swallowed.
H315	Causes skin irritation.
H317	May cause an allergic skin reaction.
H318	Causes serious eye damage.
H319	Causes serious eye irritation.
H334	May cause allergy or asthma symptoms or breathing difficulties if inhaled.
H335	May cause respiratory irritation.
H341	Suspected of causing genetic defects.
H350i	May cause cancer by inhalation.
H360F	May damage fertility.
H360FD	May damage fertility. May damage the unborn child.
H373	May cause damage to organs through prolonged or repeated exposure.
H400	Very toxic to aquatic life.
H410	Very toxic to aquatic life with long lasting effects
H411	Toxic to aquatic life with long lasting effects.

SDS Biochemicals version 2018

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product