

FICHE TECHNIQUE



A0202: ANDERSON'S RHODODENDRON INCLUDING VITAMINS

Micro Elements

mg/l

μM

CoCl₂.6H₂O

0.025

0.11

CuSO₄.5H₂O

0.025

0.10

FeNaEDTA

73.40

200.00

H₃BO₃

6.20

100.27

KI

0.30

1.81

MnSO₄.H₂O

16.90

100.00

Na₂MoO₄.2H₂O

0.25

1.03

ZnSO₄.7H₂O

8.60

29.91

Macro Elements

mg/l

mM

CaCl₂

332.02

2.99

KNO3
480.00
4.75
MgSO4
180.54
1.50
NaH2PO4
330.60
2.75
NH4NO3
400.00
5.00

Vitamins

mg/l

µM

Adenine sulphate

80.00

197.87

myo-Inositol

100.00

554.94

Thiamine HCl

0.40

1.19

Total concentration Micro and Macro elements including -vitamins: 2009.26 mg/l
Anderson W.C., Tissue culture propagation of Red and Black Raspberries, *Rubus Idaeus*
and *R. Occidentalis* Act. Hort., 112, 13 (1980).

Hazard statements

- **H315** Causes skin irritation
- **H319** Causes serious eye irritation
- **H335** May cause respiratory irritation

Precautionary statements

- **P261** Avoid breathing dust/fume/gas/mist/vapours/spray
- **P302** IF ON SKIN:
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes
- **P352** Wash with soap and water

A0202.0001/ 1 l

A0202.0010/ 10 l

A0202.0050/ 50 l