

## FICHE TECHNIQUE



### **A0528: AMMONIUM CHLORIDE**

CAS number 12125-02-9  
Molecular weight  $\text{NH}_4\text{Cl}$  = 53.5  
Assay > 99 %

#### **Storage condition**

- Store at room temperature
- Soluble in water (760 g/l / 20 °C)

#### **Hazard statements**

- **H302** Harmful if swallowed
- **H319** Causes serious eye irritation

#### **Precautionary statements**

- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P301** IF SWALLOWED:
- **P305** IF IN EYES:
- **P312** Call a POISON CENTER or doctor/physician if you feel unwell
- **P351** Rinse continuously with water for several minutes

A0528.1000/ 1 kg