

FICHE TECHNIQUE



A0532: ALUMINIUM CHLORIDE HEXAHYDRATE

CAS number 7784-13-6

Molecular weight $\text{AlCl}_3 \cdot 6\text{H}_2\text{O} = 241.4$

Assay > 98 %

Storage condition

- Soluble in water (1330 g/l / 20 °C)

Hazard statements

- **H315** Causes skin irritation
- **H319** Causes serious eye irritation

Precautionary statements

- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes

A0532.0025/ 25 g

A0532.0100/ 100 g