

## FICHE TECHNIQUE



### **A0704: L-ARGININE**

CAS number 74-79-3

Molecular weight  $C_6H_{14}N_4O_2 = 174.2$

Assay > 98.5 %

Foreign amino acids < 0.3 %

#### **Storage condition**

- Store at room temperature
- Soluble in water (150 g/l / 20 °C)

#### **Hazard statements**

- **H319** Causes serious eye irritation

#### **Precautionary statements**

- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P305** IF IN EYES:
- **P313** Get medical advice/attention
- **P337** If eye irritation persists:
- **P351** Rinse continuously with water for several minutes

A0704.0025/ 25 g

A0704.0100/ 100 g

A0704.0500/ 500 g

A0704.1000/ 1 kg