

## FICHE TECHNIQUE



### **A0941: ABSCISIC ACID (S-ABA)**

CAS number 21293-29-8

Molecular weight C<sub>15</sub>H<sub>20</sub>O<sub>4</sub> = 264.3

Assay (HPLC) > 98 %

[α]<sub>D</sub><sup>20</sup> = +425° (c= 0.052, MeOH)

#### **Storage condition**

- Store at -25 to -15 °C

#### **Hazard statements**

- **H315** Causes skin irritation
- **H319** Causes serious eye irritation
- **H335** May cause respiratory irritation

#### **Precautionary statements**

- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P302** IF ON SKIN:
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes
- **P352** Wash with soap and water

A0941.0100/ 100 mg

A0941.0250/ 250 mg

A0941.1000/ 1 g