

FICHE TECHNIQUE

A1203: AGAROSE SPI

CAS number 9012-36-6
Gel strength 1 % > 1200 g/cm²
Gel strength 1.5 % > 2500 g/cm²
Gelling temperature 34.5 - 37.5 C.
Melting temperature 86.5 89.5 C.
Sulphate < 0.2 %
Elektroendosmosis 0.09 - 0.13
Residue on ignition < 0.5 %
Loss on drying < 7 %
DNA Binding None detected
DNAse and RNAse activity None detected

Storage condition

- Store at room temperature
- Soluble in water
- Gel strength: > 2500 g/cm²

A1203.0100/ 100 g
A1203.0500/ 500 g
A1203.1000/ 1 kg