

## FICHE TECHNIQUE

### **A1335: ADENOSINE-5-TRIPHOSPHATE**

CAS number 987-65-5  
Molecular weight  $C_{10}H_{14}N_5Na_2O_{13}P_3 \cdot 3H_2O = 605.2$   
Assay (dry weight) > 96 %  
Dry weight > 90 %  
Heavy metals < 0,002 %  
White crystalline powder

#### **Storage condition**

- Store dry at 2-8 °C
- Soluble in water

A1335.0001/ 1 g

A1335.0005/ 5 g

A1335.0010/ 10 g