

## FICHE TECHNIQUE



### **A1338: AMMONIUM DIHYDROGEN PHOSPHATE**

CAS number 7722-76-1

Molecular weight (NH<sub>4</sub>)H<sub>2</sub>PO<sub>4</sub> = 115.0

Assay > 99 %

#### **Storage condition**

- Store at room temperature
- Soluble in water (370 g/l / 20 °C)

#### **Hazard statements**

- **H319** Causes serious eye irritation

#### **Precautionary statements**

- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes

A1338.1000/ 1 kg