

D(+)-Biotin

Safety Data Sheet

B0603

Version: 2.0

according to Regulation (EU) 2015/830

Supersedes 08/08/2011

Date of issue: 20/04/2011 Revision date: 26/09/2017

SECTION 1: Identification of the substance/mixture and of the company/undertaking

Product identifier

Product form : Substance Trade name : D(+)-Biotin EC-No. : 200-399-3 CAS-No. : 58-85-5 Product code : B0603

Formula : C10H16N2O3S

1.2. Relevant identified uses of the substance or mixture and uses advised against

1.2.1. Relevant identified uses

Main use category : Professional use

Industrial/Professional use spec : For professional use only. Duchefa Biochemie B.V. products are intended

only for "in vitro laboratory" research purposes.

1.2.2. Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

Supplier

Duchefa Biochemie B.V. A. Hofmanweg 71 2031 BH Haarlem - The Netherlands T +31(0)23-5319093 - F +31(0)23-5318027 info@duchefa.nl

1.4. **Emergency telephone number**

Emergency number : Supplier contact information:

+31(0)23-5319093 (M-F 09:00-17:00) +31(0)6-30109355 (outside office hours)

SDS Ref.: B0603

Organisation/Company	Address	Comment
World Health Organization world	http://apps.who.int/poisoncentres/	Consult website for a local poison centre
directory of poison centres		

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Not classified

Adverse physicochemical, human health and environmental effects

No additional information available

2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

No labelling applicable

2.3. Other hazards

No additional information available

SECTION 3: Composition/information on ingredients

3.1. Substances

 Name
 : D(+)-Biotin

 CAS-No.
 : 58-85-5

 EC-No.
 : 200-399-3

Name	Product identifier	%
Biotin	(CAS-No.) 58-85-5	> 97,5
	(EC-No.) 200-399-3	

Full text of H-statements: see section 16

3.2. Mixtures

Not applicable

SECTION 4: First aid measures

4.1. Description of first aid measures

First-aid measures general : Seek medical attention if ill effect or irritation develops.

First-aid measures after inhalation : Assure fresh air breathing.

First-aid measures after skin contact : Wash skin with mild soap and water.

First-aid measures after eye contact : Rinse with water. First-aid measures after ingestion : Rinse mouth.

4.2. Most important symptoms and effects, both acute and delayed

No additional information available

4.3. Indication of any immediate medical attention and special treatment needed

No additional information available

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media : Water spray. Alcohol resistant foam. Dry chemical powder. Carbon dioxide

(CO2).

5.2. Special hazards arising from the substance or mixture

Hazardous decomposition products in : When heated to decomposition, emits toxic fumes. COx. NOx. SOx.

case of fire

5.3. Advice for firefighters

Protection during firefighting : Wear proper protective equipment.

Other information : Prevent fire fighting water from entering the environment.

according to Regulation (EU) 2015/830

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

General measures : Wear suitable protective clothing.

6.1.1. For non-emergency personnel

No additional information available

6.1.2. For emergency responders

No additional information available

6.2. Environmental precautions

Prevent entry to sewers and public waters.

6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Dispose in a safe manner in accordance with local/national regulations.

Sweep up dry powder and dispose properly.

6.4. Reference to other sections

No additional information available

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions for safe handling : Local exhaust is recommended where dust may occur.

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Store in dry, cool, well-ventilated area. Keep container tightly closed and

dry.

Storage temperature : 2 - 8 °C

7.3. Specific end use(s)

For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

No additional information available

8.2. Exposure controls

Hand protection:

Туре	Material	Permeation	Thickness (mm)	Standard
Gloves	Nitrile rubber (NBR)	6 (> 480 minutes)	0,11	EN 374

Eye protection:

Chemical goggles or safety glasses

Skin and body protection:

Wear suitable protective clothing

Respiratory protection:

Where excessive dust may result, wear approved mask. Type P1 (EN 143)

SDS Ref.: B0603

according to Regulation (EU) 2015/830

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state : Solid

Appearance : Powder. Crystalline powder.

Molecular mass : 244,31 g/mol Colour : White.
Odour : Odourless.

Odour threshold : No data available pH : No data available Relative evaporation rate : No data available

(butylacetate=1)

Relative density

: 229 - 232 °C Melting point Freezing point : No data available **Boiling point** : No data available Flash point : No data available : No data available Auto-ignition temperature : No data available Decomposition temperature : No data available Flammability (solid, gas) Vapour pressure : No data available Relative vapour density at 20 °C : No data available

Solubility : Poorly soluble in water. Acidic aqueous solution. Alkaline mixture.

: No data available

Water: 0,2 g/l (at 90 °C)

Log Pow : No data available
Viscosity, kinematic : No data available
Viscosity, dynamic : No data available
Explosive properties : No data available
Oxidising properties : No data available
Explosive limits : No data available

9.2. Other information

No additional information available

SECTION 10: Stability and reactivity

10.1. Reactivity

Stable under normal conditions.

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

No dangerous reactions known under normal conditions of use.

10.4. Conditions to avoid

Moisture. Heat.

10.5. Incompatible materials

No additional information available

SDS Ref.: B0603

D(+)-Biotin

Safety Data Sheet

according to Regulation (EU) 2015/830

10.6. Hazardous decomposition products

May liberate toxic gases. COx. NOx. SOx.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity : Not classified

D(+)-Biotin (58-85-5)	
LD50 oral rat	> 2000 mg/kg
Skin corrosion/irritation	: Not classified
Serious eye damage/irritation	: Not classified
Respiratory or skin sensitisation	: Not classified

Germ cell mutagenicity : Not classified Carcinogenicity : Not classified

Reproductive toxicity : Not classified STOT-single exposure : Not classified STOT-repeated exposure : Not classified

Aspiration hazard : Not classified

SECTION 12: Ecological information

12.1. Toxicity

No additional information available

12.2. Persistence and degradability

D(+)-Biotin (58-85-5)	
Persistence and degradability	Minimally biodegradable.

12.3. Bioaccumulative potential

No additional information available

12.4. Mobility in soil

No additional information available

12.5. Results of PBT and vPvB assessment

No additional information available

12.6. Other adverse effects

Additional information : Prevent entry to sewers and public waters. Avoid release to the

environment.

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Waste treatment methods : Avoid release to the environment.

Product/Packaging disposal : Disposal must be done according to official regulations.

recommendations

SDS Ref.: B0603

according to Regulation (EU) 2015/830

SECTION 14: Transport information

In accordance with ADR / IATA / IMDG

ADR	IMDG	IATA	
14.1. UN number			
Not regulated	Not regulated	Not regulated	
14.2. UN proper shipping name			
Not regulated	Not regulated	Not regulated	
Not regulated	Not regulated	Not regulated	
14.3. Transport hazard class(es)			
Not regulated	Not regulated	Not regulated	
Not regulated	Not regulated	Not regulated	
14.4. Packing group			
Not regulated	Not regulated	Not regulated	
14.5. Environmental hazards			
Not regulated	Not regulated	Not regulated	
No supplementary information available			

14.6. Special precautions for user

. Overland transport

Not regulated

- Transport by sea

Not regulated

- Air transport

Not regulated

14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

15.1.1. EU-Regulations

No REACH Annex XVII restrictions

D(+)-Biotin is not on the REACH Candidate List

D(+)-Biotin is not on the REACH Annex XIV List

15.1.2. National regulations

Ensure all national/local regulations are observed.

Germany

VwVwS Annex reference : Water hazard class (WGK) 1, low hazard to waters (Classification according

to VwVwS, Annex 3; ID No. 5330)

12th Ordinance Implementing the Federal Immission Control Act -

: Is not subject of the 12. BlmSchV (Hazardous Incident Ordinance)

SDS Ref.: B0603

12.BImSchV

Netherlands

D(+)-Biotin

Safety Data Sheet

according to Regulation (EU) 2015/830

SZW-lijst van kankerverwekkende

stoffen

SZW-lijst van mutagene stoffen

NIET-limitatieve lijst van voor de voortplanting giftige stoffen -

Borstvoeding

NIET-limitatieve lijst van voor de voortplanting giftige stoffen -Vruchtbaarheid

NIET-limitatieve lijst van voor de voortplanting giftige stoffen -Ontwikkeling

: The substance is not listed

: The substance is not listed : The substance is not listed

: The substance is not listed

: The substance is not listed

Chemical safety assessment 15.2.

No additional information available

SECTION 16: Other information

Indication of changes:

1.4	Emergency number	Modified	
8.2	Hand protection	Modified	Specified material, thickness, et
			cetera of gloves

Abbreviations and acronyms:

ATE	Acute Toxicity Estimate
ADR	European Agreement concerning the International Carriage of Dangerous Goods by Road
BCF	Bioconcentration factor
CLP	Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008
DPD	Dangerous Preparations Directive 1999/45/EC
DSD	Dangerous Substances Directive 67/548/EEC
IATA	International Air Transport Association
IMDG	International Maritime Dangerous Goods
LC50	Median lethal concentration
LD50	Median lethal dose
LOAEL	Lowest Observed Adverse Effect Level
NOAEC	No-Observed Adverse Effect Concentration
PBT	Persistent Bioaccumulative Toxic
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006
SDS	Safety Data Sheet

Data sources

: Supplier. REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006.

SDS Ref.: B0603

SDS Biochemicals version 2017

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product