

## FICHE TECHNIQUE



### **B0904: 6-BENZYLAMINOPURINE (6-BAP)**

CAS number 1214-39-7

Molecular weight C<sub>12</sub>H<sub>11</sub>N<sub>5</sub> = 225.3

Assay > 99 %

#### **Storage condition**

- Store at room temperature
- Soluble in 1 N NaOH
- Slightly soluble in water (0.44 g/l / 15°C)

#### **Hazard statements**

- **H302** Harmful if swallowed
- **H315** Causes skin irritation
- **H319** Causes serious eye irritation
- **H335** May cause respiratory irritation

#### **Precautionary statements**

- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P301** IF SWALLOWED:
- **P312** Call a POISON CENTER or doctor/physician if you feel unwell

B0904.0001/ 1 g

B0904.0005/ 5 g

B0904.0025/ 25 g