Safety Data Sheet

1.1.

according to Regulation (EU) 2015/830 Date of issue: 07/06/2011 Revision date 19/06/2018 Supersedes 02/08/2011

Product identifier

B0932

Version: 2.0

SECTION 1: Identification of the substance/mixture and of the company/undertaking

Product form	: Substance
Trade name	: N-Benzyl-9-(tetrahydropyranyl)-adenine (BPA)
CAS-No.	: 2312-73-4
Product code	: B0932
Formula	: C17H19N5O
Product group	: Raw material

1.2. Relevant identified uses of the substance or mixture and uses advised against

1.2.1. Relevant identified uses

Main use category	: Professional use
Industrial/Professional use spec	For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes.

1.2.2. Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

Supplier

Duchefa Biochemie B.V. A. Hofmanweg 71 2031 BH Haarlem - The Netherlands T +31(0)23-5319093 - F +31(0)23-5318027 info@duchefa.nl

1.4. Emergency telephone number

Emergency number

: Supplier contact information:

+31(0)23-5319093 (M-F 09:00-17:00)

+31(0)6-30109355 (outside office hours)

Organisation/Company	Address	Comment
World Health Organization world directory of poison centres	http://apps.who.int/poisoncentres/	Consult website for a local poison centre

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Acute toxicity (oral), Category 4

H302

Full text of H statements : see section 16

Adverse physicochemical, human health and environmental effects

Safety Data Sheet

according to Regulation (EU) 2015/830

2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

Hazard pictograms (CLP)

Signal word (CLP) Hazard statements (CLP) Precautionary statements (CLP) **!**

GHS07

: Warning

: H302 - Harmful if swallowed.

 P281 - Use personal protective equipment as required.
 P301+P312 - IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell.

2.3. Other hazards

No additional information available

SECTION 3: Composition/information on ingredients

3.1. Substances

Name CAS-No. N-Benzyl-9-(tetrahydropyranyl)-adenine (BPA)2312-73-4

Name	Product identifier	%
6-Benzyl-9-(2-tetrahydropyranyl)-adenine	(CAS-No.) 2312-73-4	>= 98,5

Full text of H-statements: see section 16

3.2. Mixtures

Not applicable

SECTION 4: First aid measures

4.1. Description of first aid measures

First-aid measures general First-aid measures after inhalation First-aid measures after skin contact	 Seek medical attention if ill effect or irritation develops. Assure fresh air breathing. After contact with skin, wash immediately and thoroughly with plenty of water.
First-aid measures after eye contact First-aid measures after ingestion	 Rinse eyes with water as a precaution. Rinse mouth Never give anything by mouth to an unconscious person Get medical advice/attention.

4.2. Most important symptoms and effects, both acute and delayed

No additional information available

4.3. Indication of any immediate medical attention and special treatment needed

Safety Data Sheet

according to Regulation (EU) 2015/830

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media

Dry chemical powder, alcohol-resistant foam, carbon dioxide (CO2).

5.2. Special hazards arising from the substance or mixture

Hazardous decomposition products in : Carbon oxides (CO, CO2). case of fire

5.3. Advice for firefighters

```
Protection during firefighting
```

: Do not enter fire area without proper protective equipment, including respiratory protection.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

: Water spray

6.1.1. For non-emergency personnel

Measures in case of dust release	: Avoid dust formation. Ensure there is adequate ventilation. Do not breathe dust.
6.1.2. For emergency responders	
Protective equipment	: For further information refer to section 8: "Exposure controls/personal protection".

6.2. Environmental precautions

Prevent entry to sewers and public waters.

6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Sweep up or shovel. Minimize generation of dust. Keep in suitable, closed containers for disposal.

6.4. Reference to other sections

Concerning disposal elimination after cleaning, see section 13.

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions for safe handling : Avoid ventil	contact with skin and eyes. A ation.	Avoid dust formation.	Ensure adequate
--	--------------------------------------	-----------------------	-----------------

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions	: Keep cool.
	Keep container tightly closed and dry.
Storage temperature	: 2 - 8 °C

7.3. Specific end use(s)

For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

- Additional information : Provide local exhaust or general room ventilation
- 8.2. Exposure controls

Safety Data Sheet

according to Regulation (EU) 2015/830

Hand protection:

Туре	Material	Permeation	Thickness (mm)	Standard
Gloves	Nitrile rubber (NBR)	6 (> 480 minutes)	0,11	EN 374

Eye protection:

Safety glasses (to European standard EN 166 or equivalent)

Skin and body protection:

Wear suitable protective clothing

Respiratory protection:

Dust production: dust mask with filter type P2

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state Appearance Molecular mass Colour	:	Solid Powder. 309,37 g/mol White to off-white.
Odour Odour threshold	:	No data available No data available
pH		No data available
Relative evaporation rate (butylacetate=1)	:	No data available
Melting point	:	No data available
Freezing point	:	No data available
Boiling point	:	No data available
Flash point	:	No data available
Auto-ignition temperature	:	No data available
Decomposition temperature	:	No data available
Flammability (solid, gas)	:	No data available
Vapour pressure	:	No data available
Relative vapour density at 20 °C	:	No data available
Relative density	:	No data available
Solubility	:	Soluble in ethanol.
Log Pow	:	No data available
Viscosity, kinematic	:	No data available
Viscosity, dynamic	:	No data available
Explosive properties	:	No data available
Oxidising properties	:	No data available
Explosive limits	:	No data available

9.2. Other information

Safety Data Sheet

according to Regulation (EU) 2015/830

SECTION 10: Stability and reactivity

10.1. Reactivity

Stable under normal conditions of storage, handling and use.

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

No additional information available

10.4. Conditions to avoid

No additional information available

10.5. Incompatible materials

Strong oxidizers.

10.6. Hazardous decomposition products

When heated to decomposition, emits dangerous fumes:

- COx
- NOx.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity

: Oral: Harmful if swallowed.

N-Benzyl-9-(tetrahydropyranyl)-adenine (BPA) (2312-73-4)		
LD50 oral rat	1640 mg/kg	
Skin corrosion/irritation	: Not classified	
Serious eye damage/irritation	: Not classified	
Respiratory or skin sensitisation	: Not classified	
Germ cell mutagenicity	: Not classified	
Carcinogenicity	: Not classified	
Reproductive toxicity	: Not classified	
STOT-single exposure	: Not classified	
STOT-repeated exposure	: Not classified	
Aspiration hazard	: Not classified	
Other information	: See actual entry in RTECS for complete information: UO7510000.	

SECTION 12: Ecological information

12.1. Toxicity

N-Benzyl-9-(tetrahydropyranyl)-adenine (BPA) (2312-73-4)	
LC50 fish 1	2,4 mg/l Oncorhynchus mykiss (Rainbow trout)

12.2. Persistence and degradability

Safety Data Sheet

according to Regulation (EU) 2015/830

12.3. Bioaccumulative potential

No additional information available

12.4. Mobility in soil

No additional information available

12.5. Results of PBT and vPvB assessment

No additional information available

12.6. Other adverse effects

Additional information

: Avoid release to the environment. Prevent entry to sewers and public waters

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Waste treatment methods

: Dispose of this material and its container to hazardous or special waste collection point. Avoid release to the environment.

SECTION 14: Transport information

In accordance with ADR / IATA / IMDG

ADR	IMDG	IATA			
14.1. UN number					
Not regulated	Not regulated	Not regulated			
14.2. UN proper shipping name					
Not regulated	Not regulated	Not regulated			
Not regulated	Not regulated	Not regulated			
14.3. Transport hazard class(es)					
Not regulated	Not regulated	Not regulated			
Not regulated	Not regulated	Not regulated			
14.4. Packing group					
Not regulated	Not regulated	Not regulated			
14.5. Environmental hazards					
Not regulated	Not regulated	Not regulated			
No supplementary information available					

14.6. Special precautions for user

. Overland transport

Not regulated

- Transport by sea

Not regulated

- Air transport

Not regulated

14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

Safety Data Sheet

according to Regulation (EU) 2015/830

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture 15.1.1. EU-Regulations

No REACH Annex XVII restrictions

N-Benzyl-9-(tetrahydropyranyl)-adenine (BPA) is not on the REACH Candidate List N-Benzyl-9-(tetrahydropyranyl)-adenine (BPA) is not on the REACH Annex XIV List

15.1.2. National regulations

Ensure all national/local regulations are observed.

Germany

12th Ordinance Implementing the Federal Immission Control Act - 12.BImSchV	: Is not subject of the 12. BImSchV (Hazardous Incident Ordinance)	
Netherlands		
SZW-lijst van kankerverwekkende stoffen	: The substance is not listed	
SZW-lijst van mutagene stoffen	: The substance is not listed	
NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Borstvoeding	: The substance is not listed	
NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Vruchtbaarheid	: The substance is not listed	
NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Ontwikkeling	: The substance is not listed	

15.2. Chemical safety assessment

No additional information available

SECTION 16: Other information

Indication of changes:

1.4	Emergency number	Modified	
2	Classification according to Directive 67/548/EEC [DSD] or 1999/45/EC [DPD]	Removed	Obsolete
8.2	Hand protection	Modified	Specified material, thickness, et cetera of gloves

Abbreviations and acronyms:

ATE	Acute Toxicity Estimate	
ADR	European Agreement concerning the International Carriage of Dangerous Goods by Road	
BCF	Bioconcentration factor	
CLP	Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008	
DPD	Dangerous Preparations Directive 1999/45/EC	
DSD	Dangerous Substances Directive 67/548/EEC	
IATA	International Air Transport Association	
IMDG	International Maritime Dangerous Goods	
LC50	Median lethal concentration	
LD50	Median lethal dose	
LOAEL	Lowest Observed Adverse Effect Level	
NOAEC	No-Observed Adverse Effect Concentration	
PBT	Persistent Bioaccumulative Toxic	
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006	
SDS	Safety Data Sheet	

according to Regulation (EU) 2015/830

Data sources

: Supplier. REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006.

Full text of H- and EUH-statements:

Acute Tox. 4 (Oral)	Acute toxicity (oral), Category 4
H302	Harmful if swallowed.

SDS Biochemicals version 2018

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product