

FICHE TECHNIQUE



B1516: BIS-TRIS BUFFER GRADE

CAS number 6976-37-0

Molecular weight $C_8H_{19}NO_5 = 209.2$

Assay > 99 %

pH range 5.8 - 7.2

Storage condition

- Store dry at room temperature
- Soluble in water

Hazard statements

- **H315** Causes skin irritation
- **H319** Causes serious eye irritation
- **H335** May cause respiratory irritation

Precautionary statements

- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P302** IF ON SKIN:
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes
- **P352** Wash with soap and water

B1516.0100/ 100 g

B1516.0500/ 500 g