


## FICHE TECHNIQUE

### C0109: CARBENICILLIN DISODIUM



CAS number 4800-94-6

Molecular weight C<sub>17</sub>H<sub>16</sub>N<sub>2</sub>Na<sub>2</sub>O<sub>6</sub>S = 422.4

Assay > 90 %

Water < 5.5 %

#### Storage condition

- Store dry at 2-8 °C
- Soluble in water
- Hygroscopic, protect from moisture

#### Hazard statements

- **H317** May cause an allergic skin reaction
- **H334** May cause allergy or asthma symptoms or breathing difficulties if inhaled

#### Precautionary statements

- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P302** IF ON SKIN:
- **P304** IF INHALED:
- **P311** Call a POISON CENTER or doctor/physician
- **P313** Get medical advice/attention
- **P333** If skin irritation or a rash occurs:
- **P340** Remove victim to fresh air and keep at rest in a position comfortable for breathing
- **P352** Wash with soap and water

C0109.0005/ 5 g

C0109.0025/ 25 g