

Revision date 21/09/2017

Carbenicillin disodium

Safety Data Sheet

Date of issue: 13/10/2010

Supersedes 03/08/2011

according to Regulation (EU) 2015/830

C0109

Version: 2.0

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier	
Product form	: Substance
Trade name	: Carbenicillin disodium
EC-No.	: 225-360-8
CAS-No.	: 4800-94-6
Product code	: C0109
Formula	: C17H16N2Na2O6S
Synonyms	: Disodium α-Carboxybenzylpenicillin Fugacillin
Product group	: Raw material

1.2. Relevant identified uses of the substance or mixture and uses advised against

1.2.1. Relevant identified uses

Main use category	: Professional use
Industrial/Professional use spec	: For professional use only. Duchefa Biochemie B.V. products are intended
	only for "in vitro laboratory" research purposes.

1.2.2. Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

Supplier

Duchefa Biochemie B.V. A. Hofmanweg 71 2031 BH Haarlem - The Netherlands T +31(0)23-5319093 - F +31(0)23-5318027 info@duchefa.nl

1.4. Emergency telephone number

Emergency number

: Supplier contact information: +31(0)23-5319093 (M-F 09:00-17:00) +31(0)6-30109355 (outside office hours)

Organisation/Company	Address	Comment
World Health Organization world	http://apps.who.int/poisoncentres/	Consult website for a local poison centre
directory of poison centres		

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Respiratory sensitisation, Category 1	H334
Skin sensitisation, Category 1	H317

Full text of H statements : see section 16

Safety Data Sheet

according to Regulation (EU) 2015/830


Adverse physicochemical, human health and environmental effects

No additional information available

2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

Hazard pictograms (CLP)



Other hazards 2.3.

No additional information available

SECTION 3: Composition/information on ingredients

3.1. **Substances**

Name	: Carbenicillin disodium
CAS-No.	: 4800-94-6
EC-No.	: 225-360-8

Name	Product identifier	%
Carbenicillin disodium	(CAS-No.) 4800-94-6	>= 90
	(CAS-NO.) 4000-94-0 (EC-No.) 225-360-8	>= 90

Full text of H-statements: see section 16

3.2. **Mixtures**

Not applicable

SECTION 4: First aid measures

4.1. **Description of first aid measures**

First-aid measures after inhalation	: Get medical advice/attention Remove victim to fresh air.		
First-aid measures after skin contact	: Obtain medical attention		
	Remove affected clothing and wa and water, followed by warm wate	sh all exposed skin area with mild soa er rinse.	ıр
First-aid measures after eye contact	: Get medical advice/attention		
	Rinse immediately with plenty of v	vater.	
First-aid measures after ingestion	: Get immediate medical advice/att	ention	
	Rinse mouth		
	Do not induce vomiting.		
21/09/2017 Duchefa Biochemie B.V.	EN (English)	SDS Ref.: C0109	2/8

Safety Data Sheet

according to Regulation (EU) 2015/830

C0109

4.2. Most important symptoms and effects, both acute and delayed

Symptoms/effects

: Exposure may produce an allergic reaction This material or its emissions may induce blood disorders and/or aggravate pre-existing blood disorders.

4.3. Indication of any immediate medical attention and special treatment needed

No additional information available

SECTION 5: Firefighting me	asures
5.1. Extinguishing media Suitable extinguishing media	: Dry chemical powder. Water spray Foam
Unsuitable extinguishing media	Carbon dioxide. : Strong water jet.
5.2. Special hazards arising from	the substance or mixture
Hazardous decomposition products in case of fire	 Under fire conditions, hazardous fumes will be present: COx.
5.3. Advice for firefighters	
Precautionary measures fire	: Do not enter fire area without proper protective equipment, including respiratory protection. Wear proper protective equipment.
Firefighting instructions	: Prevent fire fighting water from entering the environment.
SECTION 6: Accidental release	ase measures
6.1. Personal precautions, protec	ctive equipment and emergency procedures
General measures	: Spill should be handled by trained cleaning personnel properly equipped with respiratory and eye protection.
6.1.1. For non-emergency personn	el
Measures in case of dust release	: Use good housekeeping practices to avoid rendering dust airborne.
6.1.2. For emergency responders No additional information available	
6.2. Environmental precautions Prevent entry to sewers and public wat	ters.
6.3. Methods and material for co	ntainment and cleaning up
Methods for cleaning up	: Dispose in a safe manner in accordance with local/national regulations. Sweep up dry powder and dispose properly.
6.4. Reference to other sections	

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions for safe handling

: Ensure good ventilation of the work station.

Safety Data Sheet

according to Regulation (EU) 2015/830

C0109

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions	: Store in dry, cool, well-ventilated area
	Protect from moisture.
	Protect from sunlight.
Storage temperature	: 2 - 8 °C

7.3. Specific end use(s)

For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

Additional information

: Provide local exhaust or general room ventilation

8.2. Exposure controls

Hand protection:

Туре	Material	Permeation	Thickness (mm)	Standard
Gloves	Nitrile rubber (NBR)	6 (> 480 minutes)	0,11	EN 374

Eye protection:

Safety glasses

Skin and body protection:

Wear suitable protective clothing

Respiratory protection:

Dust mask. P3

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state	: Solid
Appearance	: Powder.
Molecular mass	: 422,4 g/mol
Colour	: White to off-white.
Odour	: Odourless.
Odour threshold	: No data available
pH	: No data available
Relative evaporation rate (butylacetate=1)	: No data available
Melting point Freezing point Boiling point Flash point Auto-ignition temperature Decomposition temperature Flammability (solid, gas) Vapour pressure Relative vapour density at 20 °C	 No data available

Safety Data Sheet

according to Regulation (EU) 2015/830

Relative density Solubility	:No data available :Water: 50 g/l
Log Pow Viscosity, kinematic Viscosity, dynamic Explosive properties Oxidising properties Explosive limits	 No data available

9.2. Other information

No additional information available

SECTION 10: Stability and reactivity

10.1. Reactivity

No additional information available

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

No additional information available

10.4. Conditions to avoid

No additional information available

10.5. Incompatible materials

Strong oxidizing agent.

10.6. Hazardous decomposition products

When heated to decomposition, emits dangerous fumes

- COx
- NOx
- SOx.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity

: Not classified

Carbenicillin disodium (4800-94-6)			
LD50 oral rat	> 10000 mg/kg		
LD50 oral	> 12000 mg/kg mouse		
Skin corrosion/irritation	: Not classified		
Serious eye damage/irritation	: Not classified		
Respiratory or skin sensitisation	: May cause sensitization by	inhalation and skin contact.	
Germ cell mutagenicity	: Not classified		
Carcinogenicity	: Not classified		
Reproductive toxicity	: Not classified		
STOT-single exposure	: Not classified		
STOT-repeated exposure	: Not classified		

Safety Data Sheet

according to Regulation (EU) 2015/830

Aspiration hazard

: Not classified

Other information

: See actual entry in RTECS for complete information: ON9105000.

SECTION 12: Ecological information

12.1. Toxicity

No additional information available

12.2. Persistence and degradability

No additional information available

12.3. Bioaccumulative potential

No additional information available

12.4. Mobility in soil

No additional information available

12.5. Results of PBT and vPvB assessment

No additional information available

12.6. Other adverse effects

No additional information available

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Waste treatment methods

: Avoid release to the environment. Dispose in a safe manner in accordance with local/national regulations.

SECTION 14: Transport information

In accordance with ADR / IATA / IMDG

ADR	IMDG	IATA
14.1. UN number		
Not applicable	Not applicable	Not applicable
14.2. UN proper shipping name		
Not applicable	Not applicable	Not applicable
Not applicable	Not applicable	Not applicable
14.3. Transport hazard class(es)		
Not applicable	Not applicable	Not applicable
Not applicable	Not applicable	Not applicable
14.4. Packing group		
Not applicable	Not applicable	Not applicable
14.5. Environmental hazards		
Not applicable	Not applicable	Not applicable
No supplementary information available		

14.6. Special precautions for user

Overland transport

Not applicable

Safety Data Sheet

according to Regulation (EU) 2015/830

- Transport by sea

Not applicable

- Air transport

Not applicable

14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture 15.1.1. EU-Regulations

No REACH Annex XVII restrictions Carbenicillin disodium is not on the REACH Candidate List Carbenicillin disodium is not on the REACH Annex XIV List

15.1.2. National regulations

Germany

12th Ordinance Implementing the Federal Immission Control Act - 12.BImSchV	: Is not subject of the 12. BImSchV (Hazardous Incident Ordinance)
Netherlands	
SZW-lijst van kankerverwekkende stoffen	: The substance is not listed
SZW-lijst van mutagene stoffen	: The substance is not listed
NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Borstvoeding	: The substance is not listed
NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Vruchtbaarheid	: The substance is not listed
NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Ontwikkeling	: The substance is not listed
Denmark	
Recommendations Danish Regulation	: Young people below the age of 18 years are not allowed to use the product

15.2. Chemical safety assessment

No additional information available

SECTION 16: Other information

Indication of changes:

1.4	Emergency number	Modified	
2	Classification according to Directive 67/548/EEC [DSD] or 1999/45/EC [DPD]	Removed	Obsolete
8.2	Hand protection	Modified	Specified material, thickness, et cetera of gloves

Abbreviations and acronyms:

ATE	Acute Toxicity Estimate
ADR	European Agreement concerning the International Carriage of Dangerous Goods by Road
BCF	Bioconcentration factor
CLP	Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008

Safety Data Sheet

according to Regulation (EU) 2015/830

DPD	Dangerous Preparations Directive 1999/45/EC
DSD	Dangerous Substances Directive 67/548/EEC
IATA	International Air Transport Association
IMDG	International Maritime Dangerous Goods
LC50	Median lethal concentration
LD50	Median lethal dose
LOAEL	Lowest Observed Adverse Effect Level
NOAEC	No-Observed Adverse Effect Concentration
PBT	Persistent Bioaccumulative Toxic
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006
SDS	Safety Data Sheet

Data sources

: Manufacturer. REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006.

Full text of H- and EUH-statements:

Resp. Sens. 1	Respiratory sensitisation, Category 1
Skin Sens. 1	Skin sensitisation, Category 1
H317	May cause an allergic skin reaction.
H334	May cause allergy or asthma symptoms or breathing difficulties if inhaled.

SDS Biochemicals version 2018

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product