C0110

Version: 2.0



Cephalexin monohydrate

Safety Data Sheet
according to Regulation (EU) 2015/830

Date of issue: 03/12/2010 Revision

Supersedes 09/08/2011

10 Revision date: 01/11/2017

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Product form : Substance

Trade name : Cephalexin monohydrate

EC-No. : 239-773-6 CAS-No. : 15686-71-2 Product code : C0110

Formula : C16H17N3O4S · H2O

Synonyms : Cefalexin
Product group : Raw material

1.2. Relevant identified uses of the substance or mixture and uses advised against

1.2.1. Relevant identified uses

Main use category : Professional use

Industrial/Professional use spec : For professional use only. Duchefa Biochemie B.V. products are intended

only for "in vitro laboratory" research purposes.

1.2.2. Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

Supplier

Duchefa Biochemie B.V.
A. Hofmanweg 71
2031 BH Haarlem - The Netherlands
T +31(0)23-5319093 - F +31(0)23-5318027
info@duchefa.nl

1.4. Emergency telephone number

Emergency number : Supplier contact information:

+31(0)23-5319093 (M-F 09:00-17:00) +31(0)6-30109355 (outside office hours)

| Organisation/Company | Address | Comment |
|---|------------------------------------|---|
| World Health Organization world directory of poison centres | http://apps.who.int/poisoncentres/ | Consult website for a local poison centre |

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Respiratory sensitisation, Category 1 H334
Skin sensitisation, Category 1 H317

Full text of hazard classes and H-statements: see section 16

Adverse physicochemical, human health and environmental effects

No additional information available

Safety Data Sheet

according to Regulation (EU) 2015/830

2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

Hazard pictograms (CLP) :



GHS08

Signal word (CLP) : Danger

Hazard statements (CLP) : H317 - May cause an allergic skin reaction.

H334 - May cause allergy or asthma symptoms or breathing difficulties if

inhaled.

Precautionary statements (CLP) : P281 - Use personal protective equipment as required.

P304+P340 - IF INHALED: Remove person to fresh air and keep

comfortable for breathing.

P302+P352 - IF ON SKIN: Wash with plenty of soap and water.

P333+P313 - If skin irritation or rash occurs: Get medical advice/attention.

2.3. Other hazards

No additional information available

SECTION 3: Composition/information on ingredients

3.1. Substances

Name : Cephalexin monohydrate

CAS-No. : 15686-71-2 EC-No. : 239-773-6

| Name | Product identifier | % |
|------------------------|----------------------|------|
| Cephalexin monohydrate | (CAS-No.) 15686-71-2 | > 95 |
| | (EC-No.) 239-773-6 | |

Full text of H-statements: see section 16

3.2. Mixtures

Not applicable

SECTION 4: First aid measures

4.1. Description of first aid measures

First-aid measures after inhalation : Get immediate medical advice/attention. Move to fresh air.

First-aid measures after skin contact : Get immediate medical advice/attention. Flush with plenty of water.

First-aid measures after eye contact : Seek medical attention if ill effect or irritation develops. Rinse with water.

First-aid measures after ingestion : Seek medical attention if ill effect or irritation develops. Rinse mouth out with

water.

4.2. Most important symptoms and effects, both acute and delayed

No additional information available

4.3. Indication of any immediate medical attention and special treatment needed

No additional information available

Safety Data Sheet

according to Regulation (EU) 2015/830

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media : Alcohol resistant foam. Dry chemical powder. Water spray.

5.2. Special hazards arising from the substance or mixture

Hazardous decomposition products in : When heated to decomposition, emits toxic fumes. COx. NOx. SOx.

case of fire

5.3. Advice for firefighters

Firefighting instructions : Avoid (reject) fire-fighting water to enter environment.

Protection during firefighting : Wear proper protective equipment.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

6.1.1. For non-emergency personnel

Measures in case of dust release : Use good housekeeping practices to avoid rendering dust airborne.

6.1.2. For emergency responders

No additional information available

6.2. Environmental precautions

Prevent entry to sewers and public waters.

6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Sweep up dry powder and dispose properly. Dispose in a safe manner in

accordance with local/national regulations.

6.4. Reference to other sections

No additional information available

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions for safe handling : Avoid dust formation.

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Store in dry, cool, well-ventilated area.

Storage temperature : 2 - 8 °C

7.3. Specific end use(s)

For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

No additional information available

8.2. Exposure controls

Hand protection:

Safety Data Sheet

according to Regulation (EU) 2015/830

| Туре | Material | Permeation | Thickness (mm) | Standard |
|--------|----------------------|-------------------|----------------|----------|
| Gloves | Nitrile rubber (NBR) | 6 (> 480 minutes) | 0,11 | EN 374 |

Eye protection:

Safety glasses (to European standard EN 166 or equivalent)

Skin and body protection:

Wear suitable protective clothing

Respiratory protection:

Wear approved mask. Filter type P2 (EN 143)

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state : Solid
Appearance : Powder.
Molecular mass : 365,4 g/mol
Colour : White.

Odour : No data available
Odour threshold : No data available
pH : No data available
Relative evaporation rate : No data available

(butylacetate=1)

Melting point : No data available
Freezing point : No data available
Boiling point : No data available
Flash point : No data available
Auto-ignition temperature : No data available
Decomposition temperature : No data available
Flammability (solid gap) : No data available

Flammability (solid, gas)

Vapour pressure

Relative vapour density at 20 °C

Relative density

: No data available

Solubility : Water: Completely soluble

Log Pow : 0,65

Viscosity, kinematic : No data available
Viscosity, dynamic : No data available
Explosive properties : No data available
Oxidising properties : No data available
Explosive limits : No data available

9.2. Other information

No additional information available

SECTION 10: Stability and reactivity

10.1. Reactivity

Stable under normal conditions of storage, handling and use.

Safety Data Sheet

according to Regulation (EU) 2015/830

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

No additional information available

10.4. Conditions to avoid

Heat.

10.5. Incompatible materials

Strong oxidizers. Strong acids.

10.6. Hazardous decomposition products

When heated to decomposition, emits dangerous fumes. COx. NOx. SOx.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity : Not classified

| Cephalexin monohydrate (15686-71-2) | |
|-------------------------------------|---------|
| LD50 oral rat | > 20000 |

Skin corrosion/irritation : Not classified Serious eye damage/irritation : Not classified

Respiratory or skin sensitisation : May cause allergy or asthma symptoms or breathing difficulties if inhaled.

May cause an allergic skin reaction.

Germ cell mutagenicity : Not classified Carcinogenicity : Not classified

Reproductive toxicity : Not classified STOT-single exposure : Not classified STOT-repeated exposure : Not classified

Aspiration hazard : Not classified

Other information : See actual entry in RTECS for complete information: XI0350000.

SECTION 12: Ecological information

12.1. Toxicity

No additional information available

12.2. Persistence and degradability

No additional information available

12.3. Bioaccumulative potential

| Cephalexin monohydrate (15686-71-2) | |
|-------------------------------------|------|
| Log Pow | 0,65 |

12.4. Mobility in soil

No additional information available

Safety Data Sheet

according to Regulation (EU) 2015/830

12.5. Results of PBT and vPvB assessment

No additional information available

12.6. Other adverse effects

Additional information : Avoid release to the environment. Prevent entry to sewers and public waters

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Waste treatment methods : Dispose in a safe manner in accordance with local/national regulations.

Avoid release to the environment.

SECTION 14: Transport information

In accordance with ADR / IATA / IMDG

| ADR | IMDG | IATA |
|--|----------------|----------------|
| 14.1. UN number | | |
| Not applicable | Not applicable | Not applicable |
| 14.2. UN proper shipping name | | |
| Not applicable | Not applicable | Not applicable |
| Not applicable | Not applicable | Not applicable |
| | | |
| 14.3. Transport hazard class(es) | | |
| Not applicable | Not applicable | Not applicable |
| Not applicable | Not applicable | Not applicable |
| 14.4. Packing group | | |
| Not applicable | Not applicable | Not applicable |
| 14.5. Environmental hazards | | |
| Not applicable | Not applicable | Not applicable |
| No supplementary information available | | |

14.6. Special precautions for user

Overland transport

Not applicable

- Transport by sea

Not applicable

- Air transport

Not applicable

14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

15.1.1. EU-Regulations

No REACH Annex XVII restrictions

Cephalexin monohydrate is not on the REACH Candidate List

Cephalexin monohydrate is not on the REACH Annex XIV List

6/8

Safety Data Sheet

according to Regulation (EU) 2015/830

15.1.2. National regulations

Ensure all national/local regulations are observed.

Germany

12th Ordinance Implementing the Federal Immission Control Act -

: Is not subject of the 12. BlmSchV (Hazardous Incident Ordinance)

12.BImSchV

Netherlands

SZW-lijst van kankerverwekkende

stoffen

SZW-lijst van mutagene stoffen NIET-limitatieve lijst van voor de voortplanting giftige stoffen -

Borstvoeding

NIET-limitatieve lijst van voor de voortplanting giftige stoffen -

Vruchtbaarheid

NIET-limitatieve lijst van voor de voortplanting giftige stoffen -Ontwikkeling

: The substance is not listed

: The substance is not listed : The substance is not listed

: The substance is not listed

: The substance is not listed

Denmark

Recommendations Danish Regulation : Young people below the age of 18 years are not allowed to use the product Pregnant/breastfeeding women working with the product must not be in direct contact with the product

15.2. **Chemical safety assessment**

No additional information available

SECTION 16: Other information

Indication of changes:

| 1.4 | Emergency number | Modified | |
|-----|------------------|----------|-----------------------------------|
| 8.2 | Hand protection | Modified | Specified material, thickness, et |
| | | | cetera of gloves |

Abbreviations and acronyms:

| ATE | Acute Toxicity Estimate |
|-------|---|
| ADR | European Agreement concerning the International Carriage of Dangerous Goods by Road |
| BCF | Bioconcentration factor |
| CLP | Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008 |
| DPD | Dangerous Preparations Directive 1999/45/EC |
| DSD | Dangerous Substances Directive 67/548/EEC |
| IATA | International Air Transport Association |
| IMDG | International Maritime Dangerous Goods |
| LC50 | Median lethal concentration |
| LD50 | Median lethal dose |
| LOAEL | Lowest Observed Adverse Effect Level |
| NOAEC | No-Observed Adverse Effect Concentration |
| PBT | Persistent Bioaccumulative Toxic |
| REACH | Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006 |
| SDS | Safety Data Sheet |

Safety Data Sheet

according to Regulation (EU) 2015/830

Data sources

: ECHA (European Chemicals Agency). Supplier.

Full text of H- and EUH-statements:

| Acute Tox. 4 (Dermal) | Acute toxicity (dermal), Category 4 |
|-----------------------|--|
| Acute Tox. 4 (Oral) | Acute toxicity (oral), Category 4 |
| Eye Irrit. 2 | Serious eye damage/eye irritation, Category 2 |
| Resp. Sens. 1 | Respiratory sensitisation, Category 1 |
| Skin Irrit. 2 | Skin corrosion/irritation, Category 2 |
| Skin Sens. 1 | Skin sensitisation, Category 1 |
| STOT SE 3 | Specific target organ toxicity — Single exposure, Category 3, Respiratory tract irritation |
| H302 | Harmful if swallowed. |
| H312 | Harmful in contact with skin. |
| H315 | Causes skin irritation |
| H317 | May cause an allergic skin reaction. |
| H319 | Causes serious eye irritation. |
| H334 | May cause allergy or asthma symptoms or breathing difficulties if inhaled. |
| H335 | May cause respiratory irritation. |

SDS Biochemicals version 2017

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product

01/11/2017 Duchefa Biochemie B.V. EN (English) SDS Ref.: C0110 8/8