

FICHE TECHNIQUE

C0110 : CEPHALEXIN MONOHYDRATE

CAS number 15686-71-2

Molecular weight C₁₆H₁₇N₃O₄S.H₂O = 365.4

Assay > 95 %

Storage condition

- Store at 2-8 °C
- Soluble in water

Hazard statements

- **H317** May cause an allergic skin reaction
- **H334** May cause allergy or asthma symptoms or breathing difficulties if inhaled

Precautionary statements

- **P281** Use personal protective equipment as required
- **P302** IF ON SKIN:
- **P304** IF INHALED:
- **P313** Get medical advice/attention
- **P333** If skin irritation or a rash occurs:
- **P340** Remove victim to fresh air and keep at rest in a position comfortable for breathing
- **P352** Wash with soap and water

C0110.0005/ 5 g

C0110.0010/ 10 g