

## FICHE TECHNIQUE

### **C0116: CHLORTETRACYCLINE HCL**

CAS number 64-72-2  
Molecular weight  $C_{22}H_{23}ClN_2O_8.HCl = 515.3$   
Assay > 89.5 %



#### **Storage condition**

- Store at 2-8 °C
- Soluble in water
- Protect from light

#### **Hazard statements**

- **H315** Causes skin irritation
- **H319** Causes serious eye irritation
- **H335** May cause respiratory irritation
- **H361** Suspected of damaging fertility or the unborn child

#### **Precautionary statements**

- **P201** Obtain special instructions before use
- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P308** IF exposed or concerned:
- **P313** Get medical advice/attention

C0116.0025/ 25 g  
C0116.0100/ 100 g