

Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 25/11/2011 Revision date: 23/11/2017

C0116

Version: 2.0

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Product form : Substance

Trade name : Chlortetracycline hydrochloride

EC-No. : 200-591-7 CAS-No. : 64-72-2 Product code : C0116

Formula : C22H23CIN2O8 · HCI
Synonyms : Aureomycin hydrochloride

Product group : Raw material

1.2. Relevant identified uses of the substance or mixture and uses advised against

1.2.1. Relevant identified uses

Main use category : Professional use

Industrial/Professional use spec : For professional use only. Duchefa Biochemie B.V. products are intended

only for "in vitro laboratory" research purposes.

1.2.2. Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

Supplier

Duchefa Biochemie B.V.
A. Hofmanweg 71
2031 BH Haarlem - The Netherlands
T +31(0)23-5319093 - F +31(0)23-5318027
info@duchefa.nl

1.4. Emergency telephone number

Emergency number : Supplier contact information:

+31(0)23-5319093 (M-F 09:00-17:00) +31(0)6-30109355 (outside office hours)

Organisation/Company	Address	Comment
World Health Organization world directory of poison centres	http://apps.who.int/poisoncentres/	Consult website for a local poison centre

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Skin corrosion/irritation, Category 2
H315
Serious eye damage/eye irritation, Category 2
H319
Reproductive toxicity, Category 2
H361
Specific target organ toxicity — Single exposure, Category 3,
H335

Respiratory tract irritation

Full text of hazard classes and H-statements : see section 16

Safety Data Sheet

according to Regulation (EU) 2015/830

Adverse physicochemical, human health and environmental effects

No additional information available

2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

Hazard pictograms (CLP) :





GHS07

GHS08

Signal word (CLP) : Warning

Hazard statements (CLP) : H315 - Causes skin irritation

H319 - Causes serious eye irritation. H335 - May cause respiratory irritation.

H361 - Suspected of damaging fertility or the unborn child.

Precautionary statements (CLP) : P201 - Obtain special instructions before use.

P280 - Wear eye protection, face protection, protective clothing, protective

gloves.

P308+P313 - IF exposed or concerned: Get medical advice/attention. P305+P351 - IF IN EYES: Rinse cautiously with water for several minutes.

2.3. Other hazards

No additional information available

SECTION 3: Composition/information on ingredients

3.1. Substances

Name : Chlortetracycline hydrochloride

CAS-No. : 64-72-2 EC-No. : 200-591-7

Name	Product identifier	%
Chlortetracycline hydrochloride	(CAS-No.) 64-72-2	> 89,5
	(EC-No.) 200-591-7	

Full text of H-statements: see section 16

3.2. Mixtures

Not applicable

SECTION 4: First aid measures

4.1. Description of first aid measures

First-aid measures after inhalation : Move to fresh air. Get r

: Move to fresh air. Get medical advice/attention. Seek medical advice.

First-aid measures after skin contact : Wash skin with plenty of water and soap.

First-aid measures after eye contact : Rinse immediately with plenty of water. Get medical advice/attention.

First-aid measures after ingestion : Rinse mouth thoroughly with water. Seek medical attention if ill effect or

irritation develops.

4.2. Most important symptoms and effects, both acute and delayed

Symptoms/effects : Nausea, vomiting.

23/11/2017 Duchefa Biochemie B.V. EN (English) SDS Ref.: C0116 2/8

Safety Data Sheet

according to Regulation (EU) 2015/830

4.3. Indication of any immediate medical attention and special treatment needed

No additional information available

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media : Alcohol resistant foam. Dry chemical powder. Carbon dioxide (CO2). Water

spray.

5.2. Special hazards arising from the substance or mixture

Hazardous decomposition products in : Under fire conditions, hazardous fumes will be present. COx. NOx.

case of fire

5.3. Advice for firefighters

Protection during firefighting : Do not enter fire area without proper protective equipment, including

respiratory protection.

Other information : Prevent fire fighting water from entering the environment.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

6.1.1. For non-emergency personnel

Measures in case of dust release : Avoid dust formation.

6.1.2. For emergency responders

No additional information available

6.2. Environmental precautions

Avoid release to the environment.

6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Sweep up dry powder and dispose properly.

6.4. Reference to other sections

No additional information available

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions for safe handling : Handle in accordance with good industrial hygiene and safety practice.

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Store in dry, cool, well-ventilated area. Keep container tightly closed and

protected from light.

Storage temperature : 2 - 8 °C

7.3. Specific end use(s)

For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

No additional information available

Safety Data Sheet

according to Regulation (EU) 2015/830

8.2. Exposure controls

Hand protection:

Туре	Material	Permeation	Thickness (mm)	Standard
Gloves	Nitrile rubber (NBR)	6 (> 480 minutes)	0,11	EN 374

Eye protection:

Safety glasses (to European standard EN 166 or equivalent)

Skin and body protection:

Long sleeved protective clothing. Chemical resistant apron

Respiratory protection:

Wear approved mask. Filter type P3 (EN 143)

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state : Solid
Appearance : Powder.
Molecular mass : 515,3 g/mol
Colour : Yellow.
Odour : Negligible.

Odour threshold : No data available pH : 2,3 - 3,3 (10 g/l) Relative evaporation rate : No data available

(butylacetate=1)

Melting point : 210 - 215 °C Freezing point : No data available : No data available **Boiling point** : No data available Flash point : No data available Auto-ignition temperature : No data available Decomposition temperature Flammability (solid, gas) : No data available Vapour pressure : No data available Relative vapour density at 20 °C : No data available Relative density : No data available

Density : $\approx 1.5 \text{ g/cm}^3$

Solubility : Moderately soluble in water.

Log Pow : No data available
Viscosity, kinematic : No data available
Viscosity, dynamic : No data available
Explosive properties : No data available
Oxidising properties : No data available
Explosive limits : No data available

9.2. Other information

No additional information available

23/11/2017 Duchefa Biochemie B.V. EN (English) SDS Ref.: C0116 4/8

Safety Data Sheet

according to Regulation (EU) 2015/830

SECTION 10: Stability and reactivity

10.1. Reactivity

Stable under normal conditions.

10.2. Chemical stability

Stable under normal conditions of storage, handling and use.

10.3. Possibility of hazardous reactions

No additional information available

10.4. Conditions to avoid

No additional information available

10.5. Incompatible materials

Strong oxidizers.

10.6. Hazardous decomposition products

When heated to decomposition, emits dangerous fumes. COx. NOx.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity : Not classified

Chlortetracycline hydrochloride (64-72-2)	
LD50 oral rat	10300 mg/kg
LD50 oral	2314 mg/kg (mouse)

Skin corrosion/irritation : Causes skin irritation.

pH: 2,3 - 3,3 (10 g/l)

Serious eye damage/irritation : Causes serious eye irritation.

pH: 2,3 - 3,3 (10 g/l)

Respiratory or skin sensitisation : Not classified Germ cell mutagenicity : Not classified Carcinogenicity : Not classified

Reproductive toxicity : Suspected of damaging fertility or the unborn child.

STOT-single exposure : May cause respiratory irritation.

STOT-repeated exposure : Not classified

Aspiration hazard : Not classified

Other information : See actual entry in RTECS for complete information: QI7800000.

SECTION 12: Ecological information

12.1. Toxicity

No additional information available

12.2. Persistence and degradability

No additional information available

 23/11/2017 Duchefa Biochemie B.V.
 EN (English)
 SDS Ref.: C0116
 5/8

Safety Data Sheet

according to Regulation (EU) 2015/830

12.3. Bioaccumulative potential

No additional information available

12.4. Mobility in soil

No additional information available

12.5. Results of PBT and vPvB assessment

No additional information available

12.6. Other adverse effects

Additional information : Avoid release to the environment. Prevent entry to sewers and public waters

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Waste treatment methods

: Avoid release to the environment. When not empty dispose of this container at hazardous or special waste collection point. Dispose in a safe manner in accordance with local/national regulations.

SECTION 14: Transport information

In accordance with ADR / IATA / IMDG

ADR	IMDG	IATA
14.1. UN number		
Not regulated	Not regulated	Not regulated
14.2. UN proper shipping name		
Not regulated	Not regulated	Not regulated
Not regulated	Not regulated	Not regulated
14.3. Transport hazard class(es)		
Not regulated	Not regulated	Not regulated
Not regulated	Not regulated	Not regulated
14.4. Packing group		
Not regulated	Not regulated	Not regulated
14.5. Environmental hazards		
Not regulated	Not regulated	Not regulated
No supplementary information available		

14.6. Special precautions for user

. Overland transport

Not regulated

- Transport by sea

Not regulated

- Air transport

Not regulated

14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

Safety Data Sheet

according to Regulation (EU) 2015/830

SECTION 15: Regulatory information

Safety, health and environmental regulations/legislation specific for the substance or mixture 15.1.1. EU-Regulations

No REACH Annex XVII restrictions

Chlortetracycline hydrochloride is not on the REACH Candidate List Chlortetracycline hydrochloride is not on the REACH Annex XIV List

15.1.2. National regulations

Ensure all national/local regulations are observed.

Germany

12th Ordinance Implementing the Federal Immission Control Act -12.BlmSchV

: Is not subject of the 12. BlmSchV (Hazardous Incident Ordinance)

Netherlands

SZW-lijst van kankerverwekkende

stoffen

: The substance is not listed

: The substance is not listed

: The substance is not listed

SZW-lijst van mutagene stoffen NIET-limitatieve lijst van voor de voortplanting giftige stoffen -

Borstvoeding

: The substance is not listed

NIET-limitatieve lijst van voor de voortplanting giftige stoffen -

Vruchtbaarheid

NIET-limitatieve lijst van voor de voortplanting giftige stoffen -

Ontwikkeling

: The substance is not listed

Denmark

Recommendations Danish Regulation : Young people below the age of 18 years are not allowed to use the product Pregnant/breastfeeding women working with the product must not be in

direct contact with the product

The requirements from the Danish Working Environment Authorities regarding work with carcinogens must be followed during use and disposal

Chemical safety assessment

No additional information available

SECTION 16: Other information

Indication of changes:

1.4	Emergency number	Modified	
2	Classification according to Directive 67/548/EEC [DSD] or 1999/45/EC [DPD]	Removed	Obsolete
8.2	Hand protection	Modified	Specified material, thickness, et cetera of gloves

Abbreviations and acronyms:

ATE	Acute Toxicity Estimate
ADR	European Agreement concerning the International Carriage of Dangerous Goods by Road
BCF	Bioconcentration factor
CLP	Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008

Safety Data Sheet

according to Regulation (EU) 2015/830

DPD	Dangerous Preparations Directive 1999/45/EC
DSD	Dangerous Substances Directive 67/548/EEC
IATA	International Air Transport Association
IMDG	International Maritime Dangerous Goods
LC50	Median lethal concentration
LD50	Median lethal dose
LOAEL	Lowest Observed Adverse Effect Level
NOAEC	No-Observed Adverse Effect Concentration
PBT	Persistent Bioaccumulative Toxic
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006
SDS	Safety Data Sheet

Data sources

: ECHA (European Chemicals Agency). Supplier.

Full text of H- and EUH-statements:

Eye Irrit. 2	Serious eye damage/eye irritation, Category 2
Repr. 2	Reproductive toxicity, Category 2
Skin Irrit. 2	Skin corrosion/irritation, Category 2
STOT SE 3	Specific target organ toxicity — Single exposure, Category 3, Respiratory tract irritation
H315	Causes skin irritation
H319	Causes serious eye irritation.
H335	May cause respiratory irritation.
H361	Suspected of damaging fertility or the unborn child.

SDS Biochemicals version 2017

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product

23/11/2017 Duchefa Biochemie B.V. EN (English) SDS Ref.: C0116 8/8