Safety Data Sheet

Date of issue: 21/09/2010

Supersedes 21/09/2015

according to Regulation (EU) 2015/830

:**L** 2015/830 Revision date: 05/04/2017 C0117

Version: 2.0

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier	
Product form	: Substance
Substance name	: Clindamycin hydrochloride
EC-No.	: 244-398-6
CAS-No.	: 21462-39-5
Product code	: C0117
Formula	: C18H33CIN2O5S · CIH
Synonyms	: 7-Chloro-7-deoxylincomycin hydrochloride / Clindamycini hydrochloridum / Dalacin

1.2. Relevant identified uses of the substance or mixture and uses advised against

1.2.1. Relevant identified uses

Main use category	: Professional use
Industrial/Professional use spec	: For professional use only. Duchefa Biochemie B.V. products are intended
	only for "in vitro laboratory" research purposes

1.2.2. Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

Supplier

Duchefa Biochemie B.V. A. Hofmanweg 71 2031 BH Haarlem - The Netherlands T +31(0)23-5319093 - F +31(0)23-5318027 info@duchefa.nl

1.4. Emergency telephone number

Emergency number

: Supplier contact information:

+31(0)23-5319093 (M-F 09:00-17:00)

+31(0)6-30109355 (outside office hours)

Organisation/Company	Address	Comment
World Health Organization world directory of poison centres	http://apps.who.int/poisoncentres/	Consult website for a local poison centre

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]Skin corrosion/irritation, Category 2H315Serious eye damage/eye irritation, Category 2H319

Central cyc damage/cyc initation, Category 2	11010
Specific target organ toxicity — Single exposure, Category 3,	H335
Respiratory tract irritation	

Full text of hazard classes and H-statements : see section 16

Safety Data Sheet

according to Regulation (EU) 2015/830

Adverse physicochemical, human health and environmental effects

No additional information available

2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

Hazard pictograms (CLP)	GHS07
Signal word (CLP)	: Warning
Hazard statements (CLP)	 H315 - Causes skin irritation H319 - Causes serious eye irritation H335 - May cause respiratory irritation
Precautionary statements (CLP)	 P280 - Wear eye protection, face protection, protective clothing, protective gloves P302+P352 - IF ON SKIN: Wash with plenty of soap and water P305+P351 - IF IN EYES: Rinse cautiously with water for several minutes

2.3. Other hazards

No additional information available

SECTION 3: Composition/information on ingredients

3.1. Substances

Name	: Clindamycin hydrochloride		
CAS-No.	: 21462-39-5	: 21462-39-5	
EC-No.	: 244-398-6		
Name	Product identifier	%	
Clindamycin hydrochloride	e (CAS-No.) 21462-39-5	>=	

Full text of H-statements: see section 16

3.2. Mixtures

Not applicable

SECTION 4: First aid measures

4.1. Description of first aid measures

First-aid measures after inhalation First-aid measures after skin contact	 Assure fresh air breathing. If you feel unwell, seek medical advice. Wash skin with mild soap and water. Seek medical attention if ill effect or irritation develops. 		
First-aid measures after eye contact	: Rinse immediately with plenty of water. Seek medical attention if ill effect or irritation develops.		
First-aid measures after ingestion	: Rinse mouth. Seek medical attention if ill effect or irritation develops.		
4.2. Most important symptoms and effects, both acute and delayed			
Symptoms/effects	 Can affect blood cell formation. Damage to liver. Gastrointestinal complaints. 		
Symptoms/effects after skin contact	: Redness, pain.		

(EC-No.) 244-398-6

>= 84

Safety Data Sheet

according to Regulation (EU) 2015/830

Symptoms/effects after eye contact : Tears.

4.3. Indication of any immediate medical attention and special treatment needed

No additional information available

SECTION 5: Fire-fighting measures			
5.1. Extinguishing media			
Suitable extinguishing media	: Alcohol resistant foam. Carbon dioxide (CO2). Dry chemical powder. Water spray.		
Unsuitable extinguishing media	: Do not use a heavy water stream.		
5.2. Special hazards arising from the substance or mixture			
Hazardous decomposition products in case of fire	: When exposed to heat, may decompose liberating hazardous gases. COx. NOx. SOx.		
5.3. Advice for firefighters			
Precautionary measures fire	 Do not enter fire area without proper protective equipment, including respiratory protection. 		
Firefighting instructions	: Exercise caution when fighting any chemical fire. Use water spray or fog for cooling exposed containers.		

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

General measures

: Equip cleanup crew with proper protection.

6.1.1. For non-emergency personnel

No additional information available

6.1.2. For emergency responders

No additional information available

6.2. Environmental precautions

Prevent entry to sewers and public waters.

6.3. Methods and material for containment and cleaning up

Methods for cleaning up

: On land, sweep or shovel into suitable containers. Flush residue with large amounts of water.

6.4. Reference to other sections

No additional information available

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Additional hazards when processed : Handle in accordance with good industrial hygiene and safety procedures.

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions	: Store in dry, well-ventilated area. Store at room temperature. Keep container closed when not in use.
Storage temperature	: Store at room temperature.

Safety Data Sheet

according to Regulation (EU) 2015/830

7.3. Specific end use(s)

For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

: Provide local exhaust or general room ventilation

8.2. Exposure controls

Hand protection:

Additional information

Туре	Material	Permeation	Thickness (mm)	Standard
Gloves	Nitrile rubber (NBR)	6 (> 480 minutes)	0,11	EN 374

Eye protection:

Safety glasses

Skin and body protection:

If repeated skin contact or contamination of clothing is likely, protective clothing should be worn. Avoid repeated or prolonged skin contact

Respiratory protection:

Wear approved mask. filter type P2 (EN 143)

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state Appearance Molecular mass Colour Odour Odour Odour threshold pH Relative evaporation rate	 Solid Crystals. 461,5 g/mol White. Negligible. No data available 3 - 5 10% in H2O No data available
(butylacetate=1) Melting point Freezing point Boiling point Flash point Auto-ignition temperature Decomposition temperature Flammability (solid, gas) Vapour pressure Relative vapour density at 20 °C Relative density Solubility	 141 - 143 °C No data available So data available No data available No data available So data available
Log Pow Viscosity, kinematic	:1,015 at 20 °C :No data available

Safety Data Sheet

according to Regulation (EU) 2015/830

Viscosity, dynamic	: No data available
Explosive properties	: No data available
Oxidising properties	: No data available
Explosive limits	: No data available

9.2. Other information

No additional information available

SECTION 10: Stability and reactivity

10.1. Reactivity

No additional information available

10.2. Chemical stability

No additional information available

10.3. Possibility of hazardous reactions

Stable under normal conditions.

10.4. Conditions to avoid

Moisture.

10.5. Incompatible materials

Strong oxidizers. Strong acids.

10.6. Hazardous decomposition products

According to process conditions, hazardous decomposition products may be generated. COx. NOx. SOx.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity

: Not classified

Clindamycin hydrochloride (21462-39-5)		
LD50 oral rat	2193 mg/kg	
Skin corrosion/irritation	: Causes skin irritation. pH: 3 - 5 10% in H2O	
Serious eye damage/irritation	: Causes serious eye irritation.	
	pH: 3 - 5 10% in H2O	
Respiratory or skin sensitisation	: Not classified	
Germ cell mutagenicity	: Not classified	
Carcinogenicity	: Not classified	
Reproductive toxicity	: Not classified	
STOT-single exposure	: May cause respiratory irritation.	
STOT-repeated exposure	: Not classified	
Aspiration hazard	: Not classified	
Other information	: See actual entry in RTECS for complete information: . GF2275000.	

Safety Data Sheet

according to Regulation (EU) 2015/830

SECTION 12: Ecological information

12.1. Toxicity

No additional information available

12.2. Persistence and degradability

No additional information available

12.3. Bioaccumulative potential

Clindamycin hydrochloride (21462-39-5)	
Log Pow	1,015 at 20 °C

12.4. Mobility in soil

No additional information available

12.5. Results of PBT and vPvB assessment

No additional information available

12.6. Other adverse effects

Additional information

: Prevent entry to sewers and public waters. Avoid release to the environment

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Waste treatment methods

: Dispose of this material and its container at hazardous or special waste collection point. Dispose in a safe manner in accordance with local/national regulations.

SECTION 14: Transport information

In accordance with ADR / RID / IMDG / IATA / ADN

ADR	IMDG	ΙΑΤΑ	ADN	RID	
14.1. UN number					
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable	
14.2. UN proper shi	14.2. UN proper shipping name				
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable	
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable	
14.3. Transport haz	ard class(es)				
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable	
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable	
14.4. Packing group					
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable	
14.5. Environmental hazards					
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable	
No supplementary information available					

14.6. Special precautions for user

- Overland transport

Not applicable

- Transport by sea

Not applicable

Safety Data Sheet

according to Regulation (EU) 2015/830

- Air transport

Not applicable

- Inland waterway transport Not applicable

- Rail transport

Not applicable

14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

15.1.1. EU-Regulations

No REACH Annex XVII restrictions Clindamycin hydrochloride is not on the REACH Candidate List Clindamycin hydrochloride is not on the REACH Annex XIV List

15.1.2. National regulations

Ensure all national/local regulations are observed

Germany

12th Ordinance Implementing the Federal Immission Control Act - 12.BImSchV	: Is not subject of the 12. BImSchV (Hazardous Incident Ordinance)
Netherlands	
SZW-lijst van kankerverwekkende stoffen	: The substance is not listed
SZW-lijst van mutagene stoffen	: The substance is not listed
NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Borstvoeding	: The substance is not listed
NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Vruchtbaarheid	: The substance is not listed
NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Ontwikkeling	: The substance is not listed
Denmark	

Recommendations Danish Regulation : Young people below the age of 18 years are not allowed to use the product

15.2. Chemical safety assessment

No additional information available

SECTION 16: Other information

Indication of changes:

2.1	Classification according to Directive 67/548/EEC [DSD] or 1999/45/EC [DPD]	Removed	Obsolete
7.2	Storage temperature	Modified	Changed storage conditions from 2

Safety Data Sheet

according to Regulation (EU) 2015/830

			- 8 ° C to room temperature.
8.2	Hand protection	Modified	Specified material, thickness, et
			cetera of gloves

Abbreviations and acronyms:

European Agreement concerning the International Carriage of Dangerous Goods by Road		
Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008		
Dangerous Preparations Directive 1999/45/EC		
Dangerous Substances Directive 67/548/EEC		
Median lethal dose		
Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006		
Safety Data Sheet		
European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways		
International Air Transport Association		
International Maritime Dangerous Goods		
Persistent Bioaccumulative Toxic		

Data sources

: Supplier.

Full text of H- and EUH-statements:

Eye Irrit. 2	Serious eye damage/eye irritation, Category 2
Skin Irrit. 2	Skin corrosion/irritation, Category 2
STOT SE 3	Specific target organ toxicity — Single exposure, Category 3, Respiratory tract irritation
H315	Causes skin irritation
H319	Causes serious eye irritation
H335	May cause respiratory irritation

SDS Biochemicals version 2017

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product