

FICHE TECHNIQUE

C0118: COLISTIN SULPHATE



CAS number 1264-72-8
Potency > 19.000 units / mg

Storage condition

- Store at room temperature
- Soluble in water
- Hygroscopic

Hazard statements

- **H301** Toxic if swallowed

Precautionary statements

- **P261** Avoid breathing dust/fume/gas/mist/vapours/spray
- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P301** IF SWALLOWED:
- **P310** Immediately call a POISON CENTER or doctor/physician

C0118.0001/ 1 g
C0118.0005/ 5 g