

## FICHE TECHNIQUE



### **C0176 : CYCLOHEXIMIDE**

CAS number 66-81-9

Molecular weight C<sub>15</sub>H<sub>23</sub>NO<sub>4</sub> = 281.4

Assay > 95 %

#### **Storage condition**

- Store at 2-8 °C
- Soluble in water (< 20g/l) and ethanol

#### **Hazard statements**

- **H341** Suspected of causing genetic defects
- **H360** May damage fertility or the unborn child
- **H411** Toxic to aquatic life with long lasting effects

#### **Precautionary statements**

- **P201** Obtain special instructions before use
- **P273** Avoid release to the environment
- **P281** Use personal protective equipment as required
- **P301** IF SWALLOWED:
- **P307** IF exposed:
- **P310** Immediately call a POISON CENTER or doctor/physician
- **P313** Get medical advice/attention

C0176.0001/ 1 g

C0176.0005/ 5 g

C0176.0025/ 25 g