

Safety Data Sheet according to Regulation (EU) 2015/830 C0176

Version: 2.0

Date of issue: 19/10/2010 Revision date: 21/09/2017 Supersedes 03/08/2011

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. **Product identifier**

Product form : Substance Trade name : Cycloheximide

: cycloheximide (ISO), 4-{}{(2R)-2-[(1S,3S,5S)-3,5-dimethyl-2-oxocyclohexyl]-Chemical name

2-hydroxyethyl}}piperidine-2,6-dione

EC Index-No. : 613-140-00-8 EC-No. : 200-636-0 CAS-No. : 66-81-9 Product code : C0176 : C15H23NO4 Formula Product group : Raw material

1.2. Relevant identified uses of the substance or mixture and uses advised against

1.2.1. Relevant identified uses

Main use category : Professional use

Industrial/Professional use spec : For professional use only. Duchefa Biochemie B.V. products are intended

only for "in vitro laboratory" research purposes.

1.2.2. Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

Supplier

Duchefa Biochemie B.V. A. Hofmanweg 71 2031 BH Haarlem - The Netherlands T +31(0)23-5319093 - F +31(0)23-5318027 info@duchefa.nl

1.4. **Emergency telephone number**

: Supplier contact information: **Emergency number**

+31(0)23-5319093 (M-F 09:00-17:00) +31(0)6-30109355 (outside office hours)

SDS Ref.: C0176

Organisation/Company	Address	Comment
World Health Organization world	http://apps.who.int/poisoncentres/	Consult website for a local poison centre
directory of poison centres		

SECTION 2: Hazards identification

Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Acute toxicity (oral), Category 1 H300 H341 Germ cell mutagenicity, Category 2 Reproductive toxicity, Category 1B H360D H411 Hazardous to the aquatic environment — Chronic Hazard, Category 2

Safety Data Sheet

according to Regulation (EU) 2015/830

Full text of hazard classes and H-statements : see section 16

Adverse physicochemical, human health and environmental effects

No additional information available

2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

Hazard pictograms (CLP)







GHS06

GHS08

GHS09

Signal word (CLP) : Danger

Hazard statements (CLP) : H300 - Fatal if swallowed.

H341 - Suspected of causing genetic defects. H360D - May damage the unborn child.

H411 - Toxic to aquatic life with long lasting effects.

Precautionary statements (CLP) : P201 - Obtain special instructions before use

P281 - Use personal protective equipment as required

P301+P310 - IF SWALLOWED: immediately call a POISON CENTER or

doctor/physician

P307+P313 - IF exposed: Get medical advice/attention.

P273 - Avoid release to the environment

2.3. Other hazards

No additional information available

SECTION 3: Composition/information on ingredients

3.1. Substances

 Name
 : Cycloheximide

 CAS-No.
 : 66-81-9

 EC-No.
 : 200-636-0

 EC Index-No.
 : 613-140-00-8

Name	Product identifier	%
cycloheximide (ISO), 4-{}{(2R)-2- [(1S,3S,5S)-3,5-dimethyl-2-oxocyclohexyl]- 2-hydroxyethyl}}piperidine-2,6-dione	(CAS-No.) 66-81-9 (EC-No.) 200-636-0 (EC Index-No.) 613-140-00-8	>= 95

Full text of H-statements: see section 16

3.2. Mixtures

Not applicable

SECTION 4: First aid measures

4.1. Description of first aid measures

First-aid measures general : In case of doubt or persistent symptoms, consult always a physician.

First-aid measures after inhalation : Remove victim to fresh air. If not breathing, give artificial respiration. Get

medical advice/attention.

First-aid measures after skin contact : Take off contaminated clothing. Wash off with soap and plenty of water.

21/09/2017 Duchefa Biochemie B.V. EN (English) SDS Ref.: C0176 2/9

Safety Data Sheet

according to Regulation (EU) 2015/830

First-aid measures after eye contact : Rinse immediately with plenty of water. If eye irritation persists: Get medical

advice/attention.

First-aid measures after ingestion : Call a physician immediately. Do not induce vomiting.

4.2. Most important symptoms and effects, both acute and delayed

Symptoms/effects : Fatal if swallowed. May damage the unborn child. Damage to kidneys and

liver.

4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media : Dry powder. Carbon dioxide. Water spray.

5.2. Special hazards arising from the substance or mixture

Hazardous decomposition products in : Under fire conditions, hazardous fumes will be present.

case of fire

5.3. Advice for firefighters

Precautionary measures fire : Do not enter fire area without proper protective equipment, including

respiratory protection. Wear proper protective equipment.

Firefighting instructions : Prevent fire fighting water from entering the environment.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

6.1.1. For non-emergency personnel

Emergency procedures : Do not breathe dust.

Measures in case of dust release : Avoid dust formation.

6.1.2. For emergency responders

Protective equipment : For further information refer to section 8: "Exposure controls/personal

protection".

6.2. Environmental precautions

Notify authorities if product enters sewers or public waters. Prevent entry to sewers and public waters.

6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Sweep up dry powder and dispose properly. Minimize generation of dust.

6.4. Reference to other sections

For further information refer to section 8. Concerning disposal elimination after cleaning, see section 13.

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions for safe handling : Avoid contact with skin, eyes and clothing.

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Store in dry, well-ventilated area. Provide adequate ventilation to minimize

dust concentrations.

Storage temperature : 2 - 8 °C

Safety Data Sheet

according to Regulation (EU) 2015/830

7.3. Specific end use(s)

For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

No additional information available

8.2. Exposure controls

Appropriate engineering controls:

Ensure good ventilation of the work station.

Hand protection:

Туре	Material	Permeation	Thickness (mm)	Standard
Gloves	Polyvinylchloride (PVC)	6 (> 480 minutes)	0,35	EN 374-1, EN 374- 2, EN 374-3, EN 420

Eye protection:

Protective goggles

Skin and body protection:

Long sleeved protective clothing. Chemical resistant apron

Respiratory protection:

Wear approved mask. filter type P3 (EN 143)

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state : Solid
Appearance : Powder.
Molecular mass : 281,4 g/mol
Colour : White.

Odour : No data available
Odour threshold : No data available
pH : No data available
Relative evaporation rate : No data available

(butylacetate=1)

: 107 - 120 °C Melting point Freezing point : No data available **Boiling** point : No data available Flash point : No data available : No data available Auto-ignition temperature Decomposition temperature : No data available Flammability (solid, gas) : No data available Vapour pressure : No data available Relative vapour density at 20 °C : No data available : No data available Relative density

 21/09/2017 Duchefa Biochemie B.V.
 EN (English)
 SDS Ref.: C0176
 4/9

Safety Data Sheet

according to Regulation (EU) 2015/830

Solubility : Water: < 20 g/l

Log Pow : 0,55

Viscosity, kinematic : No data available
Viscosity, dynamic : No data available
Explosive properties : No data available
Oxidising properties : No data available
Explosive limits : No data available

9.2. Other information

No additional information available

SECTION 10: Stability and reactivity

10.1. Reactivity

No additional information available

10.2. Chemical stability

Stable under recommended handling and storage conditions (see section 7).

10.3. Possibility of hazardous reactions

No dangerous reactions known under normal conditions of use.

10.4. Conditions to avoid

No additional information available

10.5. Incompatible materials

Oxidizing agents and acids. Strong bases.

10.6. Hazardous decomposition products

When heated to decomposition, emits dangerous fumes.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity : Oral: Fatal if swallowed.

Cycloheximide (66-81-9)	
LD50 oral rat	2 mg/kg
Skin corrosion/irritation	: Not classified
Serious eye damage/irritation	: Not classified
Respiratory or skin sensitisation	: Not classified
Germ cell mutagenicity	: Suspected of causing genetic defects.
Carcinogenicity	: Not classified
Reproductive toxicity	: May damage the unborn child.
STOT-single exposure	: Not classified
STOT-repeated exposure	: Not classified
Aspiration hazard	: Not classified

SDS Ref.: C0176

Safety Data Sheet

according to Regulation (EU) 2015/830

SECTION 12: Ecological information

12.1. Toxicity

No additional information available

12.2. Persistence and degradability

No additional information available

12.3. Bioaccumulative potential

Cycloheximide (66-81-9)	
Log Pow	0,55

12.4. Mobility in soil

Cycloheximide (66-81-9)	
Ecology - soil	Prevent soil and water pollution.

12.5. Results of PBT and vPvB assessment

No additional information available

12.6. Other adverse effects

No additional information available

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Product/Packaging disposal

recommendations

: Dispose in a safe manner in accordance with local/national regulations.

Ecology - waste materials : Avoid release to the environment.

SECTION 14: Transport information

In accordance with IMDG / IATA / ADR

ADR	IMDG	IATA
14.1. UN number		
2811	2811	2811
14.2. UN proper shipping name		
TOXIC SOLID, ORGANIC, N.O.S.	TOXIC SOLID, ORGANIC, N.O.S.	Toxic solid, organic, n.o.s.
Transport document description		
UN 2811 TOXIC SOLID, ORGANIC, N.O.S., 6.1, I, (C/E), ENVIRONMENTALLY HAZARDOUS	UN 2811 TOXIC SOLID, ORGANIC, N.O.S., 6.1, I, MARINE POLLUTANT/ENVIRONMENTALLY HAZARDOUS	UN 2811 Toxic solid, organic, n.o.s., 6.1, I, ENVIRONMENTALLY HAZARDOUS
14.3. Transport hazard class(es)		
6.1	6.1	6.1
\$ P.		6
14.4. Packing group		
1	1	1
14.5. Environmental hazards		
Dangerous for the environment : Yes	Dangerous for the environment : Yes	Dangerous for the environment : Yes

Safety Data Sheet

according to Regulation (EU) 2015/830

ADR	IMDG	IATA
	Marine pollutant : Yes	
No supplementary information available		

14.6. Special precautions for user

Overland transport

Classification code (ADR) : T2

Special provisions (ADR) : 274, 614

Limited quantities (ADR) : 0 Excepted quantities (ADR) : E5 Transport category (ADR) : 1

Special provisions for carriage -

: S9, S14

: 66

Operation (ADR)

Hazard identification number (Kemler

No.)

Orange plates

66

Tunnel restriction code (ADR) : C/E EAC code : 2X

- Transport by sea

Special provisions (IMDG) : 274 Limited quantities (IMDG) : 0 Excepted quantities (IMDG) : E5 : F-A EmS-No. (Fire) EmS-No. (Spillage) : S-A Stowage category (IMDG) : B

Properties and observations (IMDG) : Toxic if swallowed, by skin contact or by inhalation.

MFAG-No : 154

- Air transport

PCA Excepted quantities (IATA) : E5

PCA Limited quantities (IATA) : Forbidden PCA limited quantity max net quantity : Forbidden

(IATA)

PCA packing instructions (IATA) : 666 PCA max net quantity (IATA) : 5kg CAO packing instructions (IATA) : 673 CAO max net quantity (IATA) : 50kg Special provisions (IATA) : A3, A5 ERG code (IATA) : 6L

Transport in bulk according to Annex II of Marpol and the IBC Code 14.7.

Not applicable

SECTION 15: Regulatory information

Safety, health and environmental regulations/legislation specific for the substance or mixture 15.1.1. EU-Regulations

SDS Ref.: C0176

Safety Data Sheet

according to Regulation (EU) 2015/830

No REACH Annex XVII restrictions

Cycloheximide is not on the REACH Candidate List Cycloheximide is not on the REACH Annex XIV List

15.1.2. National regulations

Germany

VwVwS Annex reference : Water hazard class (WGK) 3, severe hazard to waters (Classification

according to VwVwS, Annex 1 or 2; ID No. 890)

: Is not subject of the 12. BlmSchV (Hazardous Incident Ordinance)

12th Ordinance Implementing the Federal Immission Control Act -

12.BlmSchV

Netherlands

: The substance is not listed

SZW-lijst van kankerverwekkende

stoffen

SZW-lijst van mutagene stoffen NIET-limitatieve lijst van voor de

voortplanting giftige stoffen -

Borstvoeding

: The substance is not listed : The substance is not listed

NIET-limitatieve lijst van voor de voortplanting giftige stoffen -

Vruchtbaarheid

: The substance is not listed

NIET-limitatieve lijst van voor de voortplanting giftige stoffen -

Ontwikkeling

: Cycloheximide is listed

Denmark

Recommendations Danish Regulation : Young people below the age of 18 years are not allowed to use the product Pregnant/breastfeeding women working with the product must not be in

direct contact with the product

The requirements from the Danish Working Environment Authorities regarding work with carcinogens must be followed during use and disposal

SDS Ref.: C0176

15.2. **Chemical safety assessment**

No chemical safety assessment has been carried out

SECTION 16: Other information

Indication of changes:

1.4	Emergency number	Modified	
2	Classification according to Directive 67/548/EEC [DSD] or 1999/45/EC [DPD]	Removed	Obsolete

Abbreviations and acronyms:

	5 dilla dicioni,
ADN	European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways
ADR	European Agreement concerning the International Carriage of Dangerous Goods by Road
BCF	Bioconcentration factor
CLP	Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008
DPD	Dangerous Preparations Directive 1999/45/EC
DSD	Dangerous Substances Directive 67/548/EEC
IATA	International Air Transport Association
IMDG	International Maritime Dangerous Goods

Safety Data Sheet

according to Regulation (EU) 2015/830

LC50	Median lethal concentration
LD50	Median lethal dose
LOAEL	Lowest Observed Adverse Effect Level
NOAEC	No-Observed Adverse Effect Concentration
PBT	Persistent Bioaccumulative Toxic
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006
SDS	Safety Data Sheet
ATE	Acute Toxicity Estimate

Data sources

: Manufacturer. REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006.

Full text of H- and EUH-statements:

Acute Tox. 1 (Oral)	Acute toxicity (oral), Category 1
Acute Tox. 2 (Oral)	Acute toxicity (oral), Category 2
Aquatic Chronic 2	Hazardous to the aquatic environment — Chronic Hazard, Category 2
Muta. 2	Germ cell mutagenicity, Category 2
Repr. 1B	Reproductive toxicity, Category 1B
H300	Fatal if swallowed.
H341	Suspected of causing genetic defects.
H360D	May damage the unborn child.
H411	Toxic to aquatic life with long lasting effects.

SDS Biochemicals version 2017

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product

 21/09/2017 Duchefa Biochemie B.V.
 EN (English)
 SDS Ref.: C0176
 9/9