

FICHE TECHNIQUE



C0248 : CHÉE & POOL MEDIUM

Micro Elements

mg/l

μM

CoCl₂.6H₂O

0.025

0.11

CuSO₄.5H₂O

0.025

0.10

FeNaEDTA

36.70

100.00

H₃BO₃

6.20

100.27

MnSO₄.H₂O

0.85

5.00

Na₂MoO₄.2H₂O

0.25

1.03

ZnSO₄.7H₂O

8.60

29.91

Macro Elements

mg/l

mM

Ca(NO₃)₂

492.30

2.99

KH₂PO₄

170.00

1.25

KNO₃

1900.00

18.79

MgSO₄

180.54

1.50

NH₄NO₃

1650.00

20.61

Chée, R., and Pool, R.M., Improved Inorganic Media Constituents for In Vitro Shoot Multiplication of Vitis, Scientia Horticulturae, 32 (1987) 85-95.

Hazard statements

- **H315** Causes skin irritation
- **H319** Causes serious eye irritation
- **H335** May cause respiratory irritation

Precautionary statements

- **P261** Avoid breathing dust/fume/gas/mist/vapours/spray
- **P302** IF ON SKIN:
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes
- **P352** Wash with soap and water

C0248.0010/ 10 l