

Safety Data Sheet

C0401

Version: 1.0

according to Regulation (EU) 2015/830 Date of issue: 09/05/2017

# SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### 1.1. Product identifier

Product form : Mixture

Trade name : CHU vitamin mixture

Product code : C0401

Product group : Raw material

#### 1.2. Relevant identified uses of the substance or mixture and uses advised against

#### 1.2.1. Relevant identified uses

Main use category : Professional use

Industrial/Professional use spec : For professional use only. Duchefa Biochemie B.V. products are intended

only for "in vitro laboratory" research purposes.

### 1.2.2. Uses advised against

No additional information available

#### 1.3. Details of the supplier of the safety data sheet

#### Manufacturer

Duchefa Biochemie B.V.
A. Hofmanweg 71
2031 BH Haarlem - The Netherlands
T +31(0)23-5319093 - F +31(0)23-5318027
info@duchefa.nl

### 1.4. Emergency telephone number

Emergency number : Supplier contact information:

+31(0)23-5319093 (M-F 09:00-17:00) +31(0)6-30109355 (outside office hours)

Organisation/Company	Address	Comment
World Health Organization world	http://apps.who.int/poisoncentres/	Consult website for a local poison centre
directory of poison centres		

# **SECTION 2: Hazards identification**

#### 2.1. Classification of the substance or mixture

# Classification according to Regulation (EC) No. 1272/2008 [CLP]

Not classified

#### Adverse physicochemical, human health and environmental effects

No additional information available

### 2.2. Label elements

### Labelling according to Regulation (EC) No. 1272/2008 [CLP]

No labelling applicable

# 2.3. Other hazards

### No additional information available

Duchefa Biochemie B.V. EN (English) SDS Ref.: C0401 1/7

Safety Data Sheet

according to Regulation (EU) 2015/830

# SECTION 3: Composition/information on ingredients

#### 3.1. Substances

Not applicable

#### 3.2. Mixtures

Name	Product identifier	%	Classification according to Regulation (EC) No. 1272/2008 [CLP]
Glycine	(CAS-No.) 56-40-6 (EC-No.) 200-272-2	50	Not classified
Thiamine hydrochloride	(CAS-No.) 67-03-8 (EC-No.) 200-641-8	25	Not classified
Pyridoxine HCI (Vitamine B6)	(CAS-No.) 58-56-0 (EC-No.) 200-386-2	12,5	Eye Dam. 1, H318
Nicotinic acid	(CAS-No.) 59-67-6 (EC-No.) 200-441-0	12,5	Eye Irrit. 2, H319

Full text of H-statements: see section 16

# **SECTION 4: First aid measures**

### 4.1. Description of first aid measures

First-aid measures general : Seek medical attention if ill effect develops.

First-aid measures after inhalation : Remove victim to fresh air.

First-aid measures after skin contact : Wash skin with mild soap and water.

First-aid measures after eye contact : Rinse with water. First-aid measures after ingestion : Rinse mouth.

### 4.2. Most important symptoms and effects, both acute and delayed

Symptoms/effects after skin contact : May cause skin irritation.

### 4.3. Indication of any immediate medical attention and special treatment needed

No additional information available

# SECTION 5: Firefighting measures

#### 5.1. Extinguishing media

Suitable extinguishing media : Water spray

Carbon dioxide

Foam

Dry powder.

### 5.2. Special hazards arising from the substance or mixture

Hazardous decomposition products in : Under fire conditions, hazardous fumes will be present:

case of fire

- NOx

- COx

- SOx.

Duchefa Biochemie B.V. EN (English) SDS Ref.: C0401 2/7

Safety Data Sheet

according to Regulation (EU) 2015/830

### 5.3. Advice for firefighters

Precautionary measures fire : Wear proper protective equipment.

## **SECTION 6: Accidental release measures**

### 6.1. Personal precautions, protective equipment and emergency procedures

General measures : Avoid raising powdered materials into airborne dust.

#### 6.1.1. For non-emergency personnel

Emergency procedures : Wear suitable protective clothing.

### 6.1.2. For emergency responders

No additional information available

#### 6.2. Environmental precautions

Prevent entry to sewers and public waters.

#### 6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Sweep up dry powder and dispose properly.

#### 6.4. Reference to other sections

No additional information available

# **SECTION 7: Handling and storage**

### 7.1. Precautions for safe handling

Precautions for safe handling : Avoid dust formation.

# 7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Store at room temperature

Store in dry, well-ventilated area.

### 7.3. Specific end use(s)

For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes.

# SECTION 8: Exposure controls/personal protection

#### 8.1. Control parameters

No additional information available

#### 8.2. Exposure controls

#### Hand protection:

Туре	Material	Permeation	Thickness (mm)	Standard
Gloves	Nitrile rubber (NBR)	6 (> 480 minutes)	0,11	EN 374

#### Eye protection:

Safety glasses

### Skin and body protection:

In case of possible repeated skin contact wear protective clothing

#### Respiratory protection:

Where excessive dust may result, wear approved mask. Type P1 (EN 143)

Duchefa Biochemie B.V. EN (English) SDS Ref.: C0401 3/7

Safety Data Sheet

according to Regulation (EU) 2015/830

# SECTION 9: Physical and chemical properties

# 9.1. Information on basic physical and chemical properties

Physical state : Solid
Appearance : Powder.
Colour : White.

Odour : Weak.

Odour threshold : No data available pH : No data available Relative evaporation rate (butylacetate=1) : No data available

Melting point : No data available Freezing point : No data available : No data available **Boiling point** Flash point : No data available Auto-ignition temperature : No data available Decomposition temperature : No data available : No data available Flammability (solid, gas) : No data available Vapour pressure Relative vapour density at 20 °C : No data available Relative density : No data available

Solubility : Readily soluble in water.

Log Pow: No data availableViscosity, kinematic: No data availableViscosity, dynamic: No data availableExplosive properties: No data availableOxidising properties: No data availableExplosive limits: No data available

#### 9.2. Other information

No additional information available

# **SECTION 10: Stability and reactivity**

### 10.1. Reactivity

No additional information available

#### 10.2. Chemical stability

No additional information available

#### 10.3. Possibility of hazardous reactions

No additional information available

#### 10.4. Conditions to avoid

Moisture.

### 10.5. Incompatible materials

Strong oxidizing agent.

### 10.6. Hazardous decomposition products

When heated to decomposition, emits dangerous fumes:

- COx

Duchefa Biochemie B.V. EN (English) SDS Ref.: C0401 4/7

# Safety Data Sheet

according to Regulation (EU) 2015/830

- NOx
- SOx.

# **SECTION 11: Toxicological information**

### 11.1. Information on toxicological effects

Acute toxicity : Not classified

Skin corrosion/irritation : Not classified

Serious eye damage/irritation : Not classified

Respiratory or skin sensitisation : Not classified

Germ cell mutagenicity : Not classified

Carcinogenicity : Not classified

Reproductive toxicity : Not classified STOT-single exposure : Not classified STOT-repeated exposure : Not classified

Aspiration hazard : Not classified

# **SECTION 12: Ecological information**

#### 12.1. Toxicity

No additional information available

#### 12.2. Persistence and degradability

No additional information available

#### 12.3. Bioaccumulative potential

No additional information available

### 12.4. Mobility in soil

No additional information available

# 12.5. Results of PBT and vPvB assessment

No additional information available

# 12.6. Other adverse effects

Additional information : Prevent entry to sewers and public waters. Avoid release to the environment

# **SECTION 13: Disposal considerations**

#### 13.1. Waste treatment methods

Waste treatment methods : Dispose in a safe manner in accordance with local/national regulations.

Avoid release to the environment.

# **SECTION 14: Transport information**

In accordance with ADR / RID / IMDG / IATA / ADN

ADR	IMDG	IATA	ADN	RID
14.1. UN number				
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable

Duchefa Biochemie B.V. EN (English) SDS Ref.: C0401 5/7

# Safety Data Sheet

according to Regulation (EU) 2015/830

ADR	IMDG	IATA	ADN	RID	
14.2. UN proper sh	14.2. UN proper shipping name				
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable	
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable	
14.3. Transport has	zard class(es)		<u>'</u>		
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable	
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable	
14.4. Packing group					
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable	
14.5. Environmental hazards					
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable	
No supplementary information available					

#### 14.6. Special precautions for user

#### . Overland transport

Not applicable

### - Transport by sea

Not applicable

#### - Air transport

Not applicable

### - Inland waterway transport

Not applicable

### - Rail transport

Not applicable

### 14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

# **SECTION 15: Regulatory information**

# 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture 15.1.1. EU-Regulations

Contains no REACH substances with Annex XVII restrictions

Contains no substance on the REACH candidate list

Contains no REACH Annex XIV substances

### 15.1.2. National regulations

#### Germany

Reference to AwSV : Water hazard class (WGK) 1, low

hazard to water (Classification according to AwSV, Annex 1)

12th Ordinance Implementing the Federal Immission Control Act -

12.BImSchV

: Is not subject of the 12. BlmSchV (Hazardous Incident Ordinance)

### **Netherlands**

SZW-lijst van kankerverwekkende stoffen : None of the components are listed

SZW-lijst van mutagene stoffen : None of the components are listed

Duchefa Biochemie B.V. EN (English) SDS Ref.: C0401 6/7

### Safety Data Sheet

according to Regulation (EU) 2015/830

NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Borstvoeding

NIET-limitatieve lijst van voor de voortplanting giftige stoffen –

Vruchtbaarheid

NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Ontwikkeling

: None of the components are listed

: None of the components are listed

: None of the components are listed

### 15.2. Chemical safety assessment

No additional information available

# **SECTION 16: Other information**

### Abbreviations and acronyms:

ADR	European Agreement concerning the International Carriage of Dangerous Goods by Road
CLP	Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008
DPD	Dangerous Preparations Directive 1999/45/EC
DSD	Dangerous Substances Directive 67/548/EEC
LD50	Median lethal dose
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006
SDS	Safety Data Sheet
ADN	European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways
IATA	International Air Transport Association
IMDG	International Maritime Dangerous Goods
PBT	Persistent Bioaccumulative Toxic

Data sources

: REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006.

#### Full text of H- and EUH-statements:

Eye Dam. 1	Serious eye damage/eye irritation, Category 1
Eye Irrit. 2	Serious eye damage/eye irritation, Category 2
H318	Causes serious eye damage.
H319	Causes serious eye irritation.

#### SDS Biochemicals version 2018

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product

Duchefa Biochemie B.V. EN (English) SDS Ref.: C0401 7/7