

FICHE TECHNIQUE

C0401: CHU VITAMINE MIXTURE

Vitamins

mg/l

μM

Glycine

2.00

26.64

Thiamine HCl

1.00

2.96

Pyridoxine HCl

0.50

2.43

Nicotinic acid

0.50

4.06

Chu C.C, The N6 medium and its application to anther culture of -cereal crops, Proc. Symp. Plant Tissue Cult., Peking, 43 (1978).Chu C.C. et al., Establishment of an efficient medium for anther culture of rice through comparative experiments on the nitrogen - sources. Scienta Sinic., 18, 659 (1975).

C0401.0100/ >100 ml

C0401.0250/ >250 ml