

## FICHE TECHNIQUE



### **C0507: COBALT CHLORIDE HEXAHYDRATE**

CAS number 7791-13-1

Molecular weight  $\text{CoCl}_2 \cdot 6\text{H}_2\text{O}$  = 237.9

Assay > 97 %

#### **Storage condition**

- Store dry at room temperature
- Soluble in water (76 g/l / 20 °C)

#### **Hazard statements**

- **H302** Harmful if swallowed
- **H317** May cause an allergic skin reaction
- **H334** May cause allergy or asthma symptoms or breathing difficulties if inhaled
- **H341** Suspected of causing genetic defects
- **H350** May cause cancer
- **H360** May damage fertility or the unborn child
- **H410** Very toxic to aquatic life with long lasting effects

#### **Precautionary statements**

- **P201** Obtain special instructions before use
- **P273** Avoid release to the environment
- **P281** Use personal protective equipment as required
- **P304** IF INHALED:
- **P308** IF exposed or concerned:
- **P313** Get medical advice/attention
- **P341** If breathing is difficult

C0507.0025/ 25 g

C0507.0100/ 100 g