

Revision date: 28/09/2017

# Calcium carbonate

Safety Data Sheet

Date of issue: 28/09/2010

Supersedes 05/06/2012

according to Regulation (EU) 2015/830

C0529

Version: 2.0

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier	
Product form	: Substance
Trade name	: Calcium carbonate
EC-No.	: 207-439-9
CAS-No.	: 471-34-1
REACH registration No	: 01-2119486795-18
Product code	: C0529
Formula	: CaCO3
Synonyms	: Calcium monocarbonate / Carbonic acid calcium salt

#### 1.2. Relevant identified uses of the substance or mixture and uses advised against

### 1.2.1. Relevant identified uses

Main use category	: Professional use
Industrial/Professional use spec	: For professional use only. Duchefa Biochemie B.V. products are intended
	only for "in vitro laboratory" research purposes.

#### 1.2.2. Uses advised against

No additional information available

#### 1.3. Details of the supplier of the safety data sheet

Duchefa Biochemie B.V. A. Hofmanweg 71 2031 BH Haarlem - The Netherlands T +31(0)23-5319093 - F +31(0)23-5318027 info@duchefa.nl

#### 1.4. **Emergency telephone number**

Emergency number

: Supplier contact information: +31(0)23-5319093 (M-F 09:00-17:00) +31(0)6-30109355 (outside office hours)

Organisation/Company	Address	Comment
World Health Organization world directory of poison centres	http://apps.who.int/poisoncentres/	Consult website for a local poison centre

### **SECTION 2: Hazards identification**

#### Classification of the substance or mixture 2.1.

### Classification according to Regulation (EC) No. 1272/2008 [CLP] Not classified

### Adverse physicochemical, human health and environmental effects

No additional information available

Safety Data Sheet

according to Regulation (EU) 2015/830

### 2.2. Label elements

### Labelling according to Regulation (EC) No. 1272/2008 [CLP]

No labelling applicable

### 2.3. Other hazards

No additional information available

### SECTION 3: Composition/information on ingredients

### 3.1. Substances

EC-No.	: 207-439-9
CAS-No.	: 471-34-1
Name	: Calcium carbonate

Name	Product identifier	%
Calcium carbonate	(CAS-No.) 471-34-1 (EC-No.) 207-439-9	>= 98,5

Full text of H-statements: see section 16

#### 3.2. Mixtures

Not applicable

### SECTION 4: First aid measures

### 4.1. Description of first aid measures

First-aid measures after inhalation	: Remove person to fresh air and keep comfortable for breathing.
First-aid measures after skin contact	: Wash skin with mild soap and water.
First-aid measures after eye contact	: Rinse with plenty of water.
First-aid measures after ingestion	: Rinse mouth. If you feel unwell, seek medical advice.

### 4.2. Most important symptoms and effects, both acute and delayed

No additional information available

### 4.3. Indication of any immediate medical attention and special treatment needed

No additional information available

### SECTION 5: Firefighting measures

<b>5.1. Extinguishing media</b> Suitable extinguishing media	: Dry chemical, CO2, or water spray or regular foam.
<b>5.2.</b> Special hazards arising from Hazardous decomposition products in case of fire	<b>the substance or mixture</b> : Thermal decomposition generates : COx.
5.3. Advice for firefighters Protection during firefighting	: Use self-contained breathing apparatus and chemically protective clothing.

Safety Data Sheet

according to Regulation (EU) 2015/830

## SECTION 6: Accidental release measures

#### 6.1. Personal precautions, protective equipment and emergency procedures

General measures

: Avoid dust formation.

#### 6.1.1. For non-emergency personnel

Protective equipment

: Wear suitable protective clothing.

#### 6.1.2. For emergency responders

No additional information available

#### 6.2. Environmental precautions

Prevent entry to sewers and public waters.

#### 6.3. Methods and material for containment and cleaning up

Methods for cleaning up

: Sweep up, shovel or vacuum.

# 6.4. Reference to other sections

No additional information available

### SECTION 7: Handling and storage

#### 7.1. Precautions for safe handling

Precautions for safe handling	: Use good housekeeping practices during storage, transfer, handling, to
	avoid excessive dust accumulation.

#### 7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Store at room temperature. Store in a well-ventilated place. Keep container tightly closed.

### 7.3. Specific end use(s)

For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes.

### **SECTION 8: Exposure controls/personal protection**

### 8.1. Control parameters

No additional information available

#### 8.2. Exposure controls

#### Hand protection:

Туре	Material	Permeation	Thickness (mm)	Standard
Gloves	Nitrile rubber (NBR)	6 (> 480 minutes)	0,11	EN 374

#### Eye protection:

Safety glasses

#### Skin and body protection:

Wear suitable protective clothing

### **Respiratory protection:**

Dust production: dust mask with filter type P2

Safety Data Sheet

according to Regulation (EU) 2015/830

## **SECTION 9: Physical and chemical properties**

## 9.1. Information on basic physical and chemical properties

Physical state	: Solid
Appearance	: Powder.
Molecular mass	: 100,1 g/mol
Colour	: White.
Odour	: Odourless.
Odour threshold	: No data available
рН	: No data available
pH solution	: 9,5 (50 g/l, 20°C)
Relative evaporation rate	: No data available
(butylacetate=1)	
Melting point	: 825 °C
Freezing point	: No data available
Boiling point	: No data available
Flash point	: No data available
Auto-ignition temperature	: No data available
Decomposition temperature	: No data available
Flammability (solid, gas)	: No data available
Vapour pressure	: No data available
Relative vapour density at 20 °C	: No data available
Relative density	: No data available
Solubility	: Water: 13 mg/l
Log Pow	: No data available
Viscosity, kinematic	: No data available
Viscosity, dynamic	: No data available
Explosive properties	: No data available
Oxidising properties	: No data available
Explosive limits	: No data available

### 9.2. Other information

No additional information available

# SECTION 10: Stability and reactivity

### 10.1. Reactivity

No additional information available

### 10.2. Chemical stability

Stable under normal conditions.

## 10.3. Possibility of hazardous reactions

Reacts violently with : Reacts vigorously with strong oxidizers and acids.

### **10.4.** Conditions to avoid

Overheating.

### 10.5. Incompatible materials

No additional information available

C0529

Safety Data Sheet

according to Regulation (EU) 2015/830

### 10.6. Hazardous decomposition products

COx.

### **SECTION 11: Toxicological information**

### 11.1. Information on toxicological effects

: Not classified
6450 mg/kg
: Not classified
: See actual entry in RTECS for complete information: FF9335000.

## **SECTION 12: Ecological information**

### 12.1. Toxicity

No additional information available

### 12.2. Persistence and degradability

No additional information available

### 12.3. Bioaccumulative potential

No additional information available

### 12.4. Mobility in soil

No additional information available

### 12.5. Results of PBT and vPvB assessment

No additional information available

### 12.6. Other adverse effects

No additional information available

### SECTION 13: Disposal considerations

### 13.1. Waste treatment methods

Waste treatment methods	: Dispose in a safe manner in accordance with local/national regulations.
Ecology - waste materials	: Avoid release to the environment.

C0529

Safety Data Sheet

according to Regulation (EU) 2015/830

# **SECTION 14: Transport information**

### In accordance with ADR / IATA / IMDG

ADR	IMDG	ΙΑΤΑ		
14.1. UN number				
Not applicable	Not applicable	Not applicable		
14.2. UN proper shipping name				
Not applicable	Not applicable	Not applicable		
Not applicable	Not applicable	Not applicable		
14.3. Transport hazard class(es)				
Not applicable	Not applicable	Not applicable		
Not applicable	Not applicable	Not applicable		
14.4. Packing group				
Not applicable	Not applicable	Not applicable		
14.5. Environmental hazards				
Not applicable	Not applicable	Not applicable		
No supplementary information available				

### 14.6. Special precautions for user

### Overland transport

Not applicable

### - Transport by sea

Not applicable

### - Air transport

Not applicable

### 14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

### **SECTION 15: Regulatory information**

### 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

### 15.1.1. EU-Regulations

No REACH Annex XVII restrictions Calcium carbonate is not on the REACH Candidate List Calcium carbonate is not on the REACH Annex XIV List

### 15.1.2. National regulations

#### Germany

VwVwS Annex reference	: Water hazard class (WGK) nwg, Non-hazardous to water (Classification according to VwVwS, Annex 1 or 2; ID No. 317)
12th Ordinance Implementing the Federal Immission Control Act - 12.BImSchV	: Is not subject of the 12. BImSchV (Hazardous Incident Ordinance)
Netherlands SZW-lijst van kankerverwekkende stoffen	: The substance is not listed

### Safety Data Sheet

according to Regulation (EU) 2015/830

SZW-lijst van mutagene stoffen NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Borstvoeding	: The substance is not listed : The substance is not listed
NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Vruchtbaarheid	: The substance is not listed
NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Ontwikkeling	: The substance is not listed

### 15.2. Chemical safety assessment

No additional information available

### **SECTION 16: Other information**

Indication of changes:

1.4	Emergency number	Modified	
2	Classification according to Directive 67/548/EEC [DSD] or 1999/45/EC [DPD]	Removed	Obsolete
2.1	Classification according to Regulation (EC) No. 1272/2008 [CLP]	Removed	
2.2	Hazard statements (CLP)	Removed	Removed H318, H315, H335
8.2	Hand protection	Modified	Specified material, thickness, et cetera of gloves
8.2	Respiratory protection	Modified	P1 to P2

Abbreviations and acronyms:

Acute Toxicity Estimate
European Agreement concerning the International Carriage of Dangerous Goods by Road
Bioconcentration factor
Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008
Dangerous Preparations Directive 1999/45/EC
Dangerous Substances Directive 67/548/EEC
International Air Transport Association
International Maritime Dangerous Goods
Median lethal concentration
Median lethal dose
Lowest Observed Adverse Effect Level
No-Observed Adverse Effect Concentration
Persistent Bioaccumulative Toxic
Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006
Safety Data Sheet

Data sources

: Supplier. REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006. according to Regulation (EU) 2015/830

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product