

Safety Data Sheet C0530 according to Regulation (EU) 2015/830 Version: 2.0

Date of issue: 06/12/2011 Revision date: 29/09/2017

Supersedes 28/12/2011

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. **Product identifier**

Product form : Substance

Trade name : Calcium citrate tetrahydrate

EC-No. : 212-391-7 CAS-No. : 5785-44-4 Product code : C0530

Formula : Ca3(C6H5O7)2 · 4H2O

Product group : Raw material

Relevant identified uses of the substance or mixture and uses advised against 1.2.

1.2.1. Relevant identified uses

: Professional use Main use category

Industrial/Professional use spec : For professional use only. Duchefa Biochemie B.V. products are intended

only for "in vitro laboratory" research purposes.

1.2.2. Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

Supplier

Duchefa Biochemie B.V. A. Hofmanweg 71 2031 BH Haarlem - The Netherlands T +31(0)23-5319093 - F +31(0)23-5318027 info@duchefa.nl

1.4. **Emergency telephone number**

Emergency number : Supplier contact information:

+31(0)23-5319093 (M-F 09:00-17:00) +31(0)6-30109355 (outside office hours)

| Organisation/Company | Address | Comment |
|---------------------------------|------------------------------------|---|
| World Health Organization world | http://apps.who.int/poisoncentres/ | Consult website for a local poison centre |
| directory of poison centres | | |

SECTION 2: Hazards identification

Classification of the substance or mixture 2.1.

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Not classified

Adverse physicochemical, human health and environmental effects

No additional information available

2.2. Label elements

According to EC directives or the corresponding national regulations there is no labelling obligation for this product.

29/09/2017 Duchefa Biochemie B.V. EN (English) 1/7

Safety Data Sheet

according to Regulation (EU) 2015/830

2.3. Other hazards

Other hazards not contributing to the

classification

: None under normal conditions.

SECTION 3: Composition/information on ingredients

3.1. Substances

Name : Calcium citrate tetrahydrate

CAS-No. : 5785-44-4 EC-No. : 212-391-7

| Name | Product identifier | % |
|------------------------------|---------------------|------|
| Calcium citrate tetrahydrate | (CAS-No.) 5785-44-4 | > 98 |
| | (EC-No.) 212-391-7 | |

Full text of H-statements: see section 16

3.2. Mixtures

Not applicable

SECTION 4: First aid measures

4.1. Description of first aid measures

First-aid measures after inhalation : Assure fresh air breathing. Get medical advice and attention if you feel

unwell.

First-aid measures after skin contact : No irritant effect.

First-aid measures after eye contact : Ensure adequate flushing of eyes by separating eyelids with the fingers.

Consult an ophtalmologist if irritation persists.

First-aid measures after ingestion : Get medical advice/attention if you feel unwell.

4.2. Most important symptoms and effects, both acute and delayed

No additional information available

4.3. Indication of any immediate medical attention and special treatment needed

No additional information available

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media : Dry chemical powder, alcohol-resistant foam, carbon dioxide (CO2). Water

spray.

5.2. Special hazards arising from the substance or mixture

Hazardous decomposition products in : Carbon oxides (CO, CO2).

case of fire

5.3. Advice for firefighters

No additional information available

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

6.1.1. For non-emergency personnel

Protective equipment : Protective clothing.

29/09/2017 Duchefa Biochemie B.V. EN (English) 2/7

Safety Data Sheet

according to Regulation (EU) 2015/830

6.1.2. For emergency responders

No additional information available

6.2. Environmental precautions

Prevent entry to sewers and public waters.

6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Mechanically recover the product.

6.4. Reference to other sections

Reference to other sections (8, 13).

SECTION 7: Handling and storage

7.1. Precautions for safe handling

No additional information available

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Keep container tightly closed. Store in dry, well-ventilated area.

Storage temperature : 15 - 25 °C

7.3. Specific end use(s)

For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

No additional information available

8.2. Exposure controls

Hand protection:

| Туре | Material | Permeation | Thickness (mm) | Standard |
|--------|----------------------|-------------------|----------------|----------|
| Gloves | Nitrile rubber (NBR) | 6 (> 480 minutes) | 0,11 | EN 374 |

Eye protection:

Safety glasses (to European standard EN 166 or equivalent)

Skin and body protection:

Wear suitable protective clothing

Respiratory protection:

Dust production: dust mask with filter type P1

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state : Solid
Appearance : Powder.
Molecular mass : 570,5 g/mol

Colour

Odour : Nearly odourless.

29/09/2017 Duchefa Biochemie B.V. EN (English) 3/7

Safety Data Sheet

according to Regulation (EU) 2015/830

Odour threshold : No data available

pH : 6 (25%)

Relative evaporation rate : No data available

(butylacetate=1)

Melting point : 120 °C

Freezing point : No data available
Boiling point : No data available
Flash point : No data available

Auto-ignition temperature : 500 °C

Decomposition temperature : No data available Flammability (solid, gas) : No data available Vapour pressure : No data available Relative vapour density at 20 °C : No data available Relative density : No data available Density : 1,831 g/cm³

Solubility : Poorly soluble in water.

Water: 0,96 g/l

Log Pow : No data available
Viscosity, kinematic : No data available
Viscosity, dynamic : No data available
Explosive properties : No data available
Oxidising properties : No data available
Explosive limits : No data available

9.2. Other information

No additional information available

SECTION 10: Stability and reactivity

10.1. Reactivity

No additional information available

10.2. Chemical stability

No additional information available

10.3. Possibility of hazardous reactions

No additional information available

10.4. Conditions to avoid

No additional information available

10.5. Incompatible materials

No additional information available

10.6. Hazardous decomposition products

Carbon oxides (CO, CO2).

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity : Not classified

29/09/2017 Duchefa Biochemie B.V. EN (English) 4/7

Safety Data Sheet

according to Regulation (EU) 2015/830

Skin corrosion/irritation : Not classified

pH: 6 (25%)

Serious eye damage/irritation : Not classified

pH: 6 (25%)

Respiratory or skin sensitisation : Not classified Germ cell mutagenicity : Not classified Carcinogenicity : Not classified

Reproductive toxicity : Not classified STOT-single exposure : Not classified STOT-repeated exposure : Not classified

Aspiration hazard : Not classified

SECTION 12: Ecological information

12.1. Toxicity

No additional information available

12.2. Persistence and degradability

No additional information available

12.3. Bioaccumulative potential

No additional information available

12.4. Mobility in soil

No additional information available

12.5. Results of PBT and vPvB assessment

No additional information available

12.6. Other adverse effects

No additional information available

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Regional legislation (waste) : Disposal must be done according to official regulations.

SECTION 14: Transport information

In accordance with ADR / IATA / IMDG

| ADR | IMDG | IATA | |
|----------------------------------|-----------------|---------------|--|
| 14.1. UN number | 14.1. UN number | | |
| Not regulated | Not regulated | Not regulated | |
| 14.2. UN proper shipping name | | | |
| Not regulated | Not regulated | Not regulated | |
| Not regulated | Not regulated | Not regulated | |
| 14.3. Transport hazard class(es) | | | |
| Not regulated | Not regulated | Not regulated | |

29/09/2017 Duchefa Biochemie B.V. EN (English) 5/7

Safety Data Sheet

according to Regulation (EU) 2015/830

| ADR | IMDG | IATA |
|--|---------------|---------------|
| Not regulated | Not regulated | Not regulated |
| 14.4. Packing group | | |
| Not regulated | Not regulated | Not regulated |
| 14.5. Environmental hazards | | |
| Not regulated | Not regulated | Not regulated |
| No supplementary information available | | |

14.6. Special precautions for user

. Overland transport

Not regulated

- Transport by sea

Not regulated

- Air transport

Not regulated

Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

SECTION 15: Regulatory information

Safety, health and environmental regulations/legislation specific for the substance or mixture 15.1.1. EU-Regulations

No REACH Annex XVII restrictions

Calcium citrate tetrahydrate is not on the REACH Candidate List Calcium citrate tetrahydrate is not on the REACH Annex XIV List

15.1.2. National regulations

Germany

VwVwS Annex reference : Water hazard class (WGK) 1, low hazard to waters (Classification according

to VwVwS, Annex 3; ID No. 2205)

12th Ordinance Implementing the Federal Immission Control Act -12.BImSchV

: Is not subject of the 12. BlmSchV (Hazardous Incident Ordinance)

Netherlands

SZW-lijst van kankerverwekkende

stoffen

: The substance is not listed

SZW-lijst van mutagene stoffen

NIET-limitatieve lijst van voor de

: The substance is not listed

voortplanting giftige stoffen -

: The substance is not listed

Borstvoeding

NIET-limitatieve lijst van voor de voortplanting giftige stoffen -

: The substance is not listed

Vruchtbaarheid

NIET-limitatieve lijst van voor de voortplanting giftige stoffen -

Ontwikkeling

: The substance is not listed

29/09/2017 Duchefa Biochemie B.V. EN (English) 6/7

Safety Data Sheet

according to Regulation (EU) 2015/830

15.2. Chemical safety assessment

No chemical safety assessment has been carried out

SECTION 16: Other information

Indication of changes:

| 1.4 | Emergency telephone | Modified | Added additional information |
|-----|--|----------|--|
| 2 | Classification according to Directive 67/548/EEC [DSD] or 1999/45/EC [DPD] | Removed | Obsolete |
| 8.2 | Hand protection | Modified | Specified material, thickness, et cetera of gloves |

Abbreviations and acronyms:

| ATE | Acute Toxicity Estimate |
|-------|---|
| ADR | European Agreement concerning the International Carriage of Dangerous Goods by Road |
| BCF | Bioconcentration factor |
| CLP | Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008 |
| DPD | Dangerous Preparations Directive 1999/45/EC |
| DSD | Dangerous Substances Directive 67/548/EEC |
| IATA | International Air Transport Association |
| IMDG | International Maritime Dangerous Goods |
| LC50 | Median lethal concentration |
| LD50 | Median lethal dose |
| LOAEL | Lowest Observed Adverse Effect Level |
| NOAEC | No-Observed Adverse Effect Concentration |
| PBT | Persistent Bioaccumulative Toxic |
| REACH | Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006 |
| SDS | Safety Data Sheet |

Data sources

: Supplier. REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006.

SDS Biochemicals version 2017

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product

29/09/2017 Duchefa Biochemie B.V. EN (English) 7/7