

## FICHE TECHNIQUE



### **C0605: CHOLINE CHLORIDE**

CAS number 67-48-1

Molecular weight  $C_5H_{14}NOCl = 139.6$

White crystals

#### **Storage condition**

- Store at room temperature
- readily soluble in water
- Hygroscopic, Under inert gas

#### **Hazard statements**

- **H315** Causes skin irritation
- **H319** Causes serious eye irritation
- **H335** May cause respiratory irritation

#### **Precautionary statements**

- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes

C0605.0100/ 100 g