

FICHE TECHNIQUE



C0706: CYSTEINE HCL MONOHYDRATE

CAS number 7048-04-6

Molecular weight $C_3H_8NO_2S \cdot H_2O = 175.6$

Assay > 98,0 %

Storage condition

- Store at room temperature
- Soluble in water

Hazard statements

- **H315** Causes skin irritation
- **H319** Causes serious eye irritation
- **H335** May cause respiratory irritation

Precautionary statements

- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P302** IF ON SKIN:
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes
- **P352** Wash with soap and water

C0706.0100/ 100 g

C0706.0500/ 500 g

C0706.1000/ 1 kg