

FICHE TECHNIQUE



C0909 : P-CHLOROPHENOXYACETIC ACID (4-CPA)

CAS number 122-88-3

Molecular weight C₈H₇ClO₃ = 186.6

Assay > 97 %

Storage condition

- Store at room temperature
- Soluble in ethanol

Hazard statements

- **H302** Harmful if swallowed

Precautionary statements

- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P301** IF SWALLOWED:
- **P310** Immediately call a POISON CENTER or doctor/physician

C0909.0025/ 25 g

C0909.0100/ 100 g

