

FICHE TECHNIQUE



C1303: CITRIC ACID MONOHYDRATE

CAS number 5949-29-1

Molecular weight $C_6H_8O_7 \cdot H_2O = 210.1$

Assay > 99.5 %

Storage condition

- Store at room temperature
- Soluble in water (1630 g/l / 20 °C)

Hazard statements

- **H315** Causes skin irritation
- **H319** Causes serious eye irritation

Precautionary statements

- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes

C1303.1000/ 1 kg