

FICHE TECHNIQUE



C1397: CETRIMIDE (N-TETRADECYL-N,N,N-TRIMETHYLAMMONIUM BROMIDE)

CAS number 8044-71-1

Molecular weight C₁₇H₃₈NBr = 336.4

Assay > 96 %

Storage condition

- Store at room temperature
- Hygroscopic, soluble in water
- CETRIMIDE

Hazard statements

- **H302** Harmful if swallowed
- **H312** Harmful in contact with skin
- **H314** Causes severe skin burns and eye damage
- **H332** Harmful if inhaled

Precautionary statements

- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P302** IF ON SKIN:
- **P305** IF IN EYES:
- **P308** IF exposed or concerned:
- **P313** Get medical advice/attention
- **P351** Rinse continuously with water for several minutes
- **P352** Wash with soap and water

C1397.0050/ 50 g

C1397.0100/ 100 g

C1397.0500/ 500 g

C1397.1000/1 kg