

FICHE TECHNIQUE



C8001: CELLULASE R-10

CAS number 9012-54-8

Beige lyophilisate

Enzyme activity > 10,000 U/g

Loss on drying < 10 %

Storage condition

- Store dry at 2-8 °C
- Loss on drying < 10%
- From Trichoderma Viride

Hazard statements

- **H334** May cause allergy or asthma symptoms or breathing difficulties if inhaled

Precautionary statements

- **P261** Avoid breathing dust/fume/gas/mist/vapours/spray
- **P304** IF INHALED:
- **P341** If breathing is difficult

C8001.0001/ 1 g

C8001.0005/ 5 g

C8001.0010/ 10 g