

FICHE TECHNIQUE



D0205: DE GREEF & JACOBS MEDIUM

Micro Elements

mg/l

μ M

CoCl₂.6H₂O

0.0025

0.01

CuSO₄.5H₂O

0.0025

0.01

FeNaEDTA

36.70

100.00

H₃BO₃

10.62

171.76

KI

1.58

9.54

MnSO₄.H₂O

1.68

9.94

Na₂MoO₄.2H₂O

0.0025

0.01

ZnSO₄.7H₂O

1.06

3.69

Macro Elements

mg/l

mM

CaCl₂

226.50

2.04

KCl

600.00

8.05

KNO₃

2000.00

19.78

MgSO₄

244.33

2.03

NaH₂PO₄

250.00

2.08

(NH₄)₂SO₄

400.00

3.03

De Greef W. and Jacobs M. In vitro culture of the sugarbeet: Description of a cell line with high regeneration capacity., Plant Science Letters., 17, 55-61 (1979).

Hazard statements

- **H315** Causes skin irritation
- **H319** Causes serious eye irritation
- **H335** May cause respiratory irritation

Precautionary statements

- **P261** Avoid breathing dust/fume/gas/mist/vapours/spray
- **P302** IF ON SKIN:
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes
- **P352** Wash with soap and water

D0205.0005/ 5 l

D0205.0050/ 50 l