

FICHE TECHNIQUE



D0247: DKW/JUGLANS MEDIUM INCLUDING VITAMINS

Micro Elements

mg/l

μM

CuSO₄.5H₂O

0.25

1.00

FeNaEDTA

44.63

121.61

H₃BO₃

4.80

77.63

MnSO₄.H₂O

33.80

200.00

Na₂MoO₄.2H₂O

0.39

1.61

ZnSO₄.7H₂O

17.00

72.19

Macro Elements

mg/l

mM

CaCl₂

112.50

1.01

Ca(NO₃)₂.2H₂O

1664.64

8.30

KH₂PO₄

265.00

1.95

K₂SO₄
1559.00
8.95
MgSO₄
361.49
3.00

NH₄NO₃
1416.00
17.70

Total concentration Micro and Macro elements: 5479.50 mg/l

Vitamins

mg/l

µM

Glycine
2.00

26.64

myo-Inositol
100.00

554.94

Nicotinic acid
1.00

8.12

Thiamine HCl
2.00

5.93

Total concentration Micro and Macro elements including -vitamins: 5584.50mg/l
river, J.A., Kuniyuki, A.H. In Vitro Propagation of Paradox walnut Rootstock, Hort. Science, 19(4), August 1984.

Hazard statements

- **H315** Causes skin irritation
- **H319** Causes serious eye irritation
- **H335** May cause respiratory irritation

Precautionary statements

- **P261** Avoid breathing dust/fume/gas/mist/vapours/spray
- **P302** IF ON SKIN:
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes
- **P352** Wash with soap and water

D0247.0001/ 1 l
D0247.0005/ 5 l
D0247.0010/ 10
D0247.0025/ 25 l
D0247.0050/ 50 l