

FICHE TECHNIQUE

D0414: DKW/JUGLANS VITAMIN MIXTURE

Vitamins

mg/l

μM

myo-Inositol

100.00

554.94

Nicotinic acid

1.00

8.12

Pyridoxine HCl

1.00

4.86

Thiamine HCl

10.00

29.65

river, J.A., Kuniyuki, A.H. In Vitro Propagation of Paradox walnut Rootstock, Hort. Science, 19(4), August 1984.

D0414.0100/ >100 ml

D0414.0250/ >250 ml