

FICHE TECHNIQUE



D0911: 2,4-DICHLOROPHENOXYACETIC ACID (2,4 D)

CAS number 94-75-7

Molecular weight $C_8H_6Cl_2O_3 = 221.0$

Assay > 96 %

Off white to tan crystals

Storage condition

- Store at room temperature
- Soluble in ethanol or 1 N NaOH

Hazard statements

- **H302** Harmful if swallowed
- **H317** May cause an allergic skin reaction
- **H318** Causes serious eye damage
- **H335** May cause respiratory irritation
- **H412** Harmful to aquatic life with long lasting effects

Precautionary statements

- **P261** Avoid breathing dust/fume/gas/mist/vapours/spray
- **P273** Avoid release to the environment
- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P308** IF exposed or concerned:
- **P313** Get medical advice/attention

D0911.0100/ 100 g

D0911.0250/ 250 g