

Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 23/09/2011 Supersedes 07/10/2011 Revision date: 21/09/2017

**D0933** Version: 2.0

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### 1.1. Product identifier

Product form : Substance

Trade name : DL-Dihydrozeatin
CAS-No. : 14894-18-9
Product code : D0933
Formula : C10H15N5O

Synonyms : 6-(4-hydroxy-3-methylbutylamino)purine / (diH)Z / DHZ / DZ

Product group : Raw material

## 1.2. Relevant identified uses of the substance or mixture and uses advised against

## 1.2.1. Relevant identified uses

Main use category : Professional use

Industrial/Professional use spec : For professional use only. Duchefa Biochemie B.V. products are intended

only for "in vitro laboratory" research purposes.

#### 1.2.2. Uses advised against

No additional information available

## 1.3. Details of the supplier of the safety data sheet

#### Supplier

Duchefa Biochemie B.V.
A. Hofmanweg 71
2031 BH Haarlem - The Netherlands
T +31(0)23-5319093 - F +31(0)23-5318027
info@duchefa.nl

## 1.4. Emergency telephone number

Emergency number : Supplier contact information:

+31(0)23-5319093 (M-F 09:00-17:00) +31(0)6-30109355 (outside office hours)

SDS Ref.: D0933

Organisation/Company	Address	Comment
World Health Organization world	http://apps.who.int/poisoncentres/	Consult website for a local poison centre
directory of poison centres	,	·

## **SECTION 2: Hazards identification**

### 2.1. Classification of the substance or mixture

### Classification according to Regulation (EC) No. 1272/2008 [CLP]

Not classified

#### Adverse physicochemical, human health and environmental effects

No additional information available

1/7

Safety Data Sheet

according to Regulation (EU) 2015/830

#### 2.2. Label elements

## Labelling according to Regulation (EC) No. 1272/2008 [CLP]

No labelling applicable

#### 2.3. Other hazards

No additional information available

## **SECTION 3: Composition/information on ingredients**

#### 3.1. Substances

Name : DL-Dihydrozeatin CAS-No. : 14894-18-9

Product identifier	%
(CAS-No.) 14894-18-9	>= 98

Full text of H-statements: see section 16

#### 3.2. Mixtures

Not applicable

## **SECTION 4: First aid measures**

## 4.1. Description of first aid measures

First-aid measures general : Seek medical attention if ill effect or irritation develops.

First-aid measures after inhalation : Remove victim to fresh air.

First-aid measures after skin contact : Wash off with soap and plenty of water.

First-aid measures after eye contact : Rinse with water. First-aid measures after ingestion : Rinse mouth.

## 4.2. Most important symptoms and effects, both acute and delayed

No additional information available

#### 4.3. Indication of any immediate medical attention and special treatment needed

No additional information available

## SECTION 5: Firefighting measures

## 5.1. Extinguishing media

Suitable extinguishing media : Water spray. Dry chemical powder, alcohol-resistant foam, carbon dioxide

(CO2).

#### 5.2. Special hazards arising from the substance or mixture

Hazardous decomposition products in : Under fire conditions, hazardous fumes will be present. COx. NOx.

case of fire

## 5.3. Advice for firefighters

Precautionary measures fire : Wear proper protective equipment.

Firefighting instructions : Prevent fire fighting water from entering the environment.

Safety Data Sheet

according to Regulation (EU) 2015/830

## **SECTION 6: Accidental release measures**

### 6.1. Personal precautions, protective equipment and emergency procedures

General measures : Avoid raising powdered materials into airborne dust.

## 6.1.1. For non-emergency personnel

Emergency procedures : Wear suitable protective clothing.

## 6.1.2. For emergency responders

No additional information available

## 6.2. Environmental precautions

Prevent entry to sewers and public waters.

## 6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Sweep up dry powder and dispose properly. Keep in suitable, closed

containers for disposal.

#### 6.4. Reference to other sections

No additional information available

## SECTION 7: Handling and storage

## 7.1. Precautions for safe handling

Precautions for safe handling : Ensure good ventilation of the work station. Avoid dust formation.

## 7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Store in dry, well-ventilated area. Keep container tightly closed. Storage temperature : 2 - 8 °C in powder form (in liquid form store at -25 to -15 °C)

## 7.3. Specific end use(s)

For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes.

## SECTION 8: Exposure controls/personal protection

#### 8.1. Control parameters

No additional information available

#### 8.2. Exposure controls

#### Hand protection:

Туре	Material	Permeation	Thickness (mm)	Standard
Gloves	Nitrile rubber (NBR)	6 (> 480 minutes)	0,11	EN 374

### Eye protection:

Safety glasses

## Skin and body protection:

In case of possible repeated skin contact wear protective clothing

#### Respiratory protection:

Dust production: dust mask with filter type P1 (EN 143)

SDS Ref.: D0933

Safety Data Sheet

according to Regulation (EU) 2015/830

## SECTION 9: Physical and chemical properties

## 9.1. Information on basic physical and chemical properties

Physical state : Solid
Appearance : Powder.
Molecular mass : 221,3 g/mol

Colour : White. light yellow.
Odour : No data available
Odour threshold : No data available
pH : No data available
Relative evaporation rate : No data available

(butylacetate=1)

Melting point : 170 °C

Freezing point : No data available **Boiling point** : No data available Flash point : No data available Auto-ignition temperature : No data available Decomposition temperature : No data available : No data available Flammability (solid, gas) Vapour pressure : No data available Relative vapour density at 20 °C : No data available Relative density : No data available : Soluble in ethanol. Solubility

Log Pow : 0,466

Viscosity, kinematic : No data available
Viscosity, dynamic : No data available
Explosive properties : No data available
Oxidising properties : No data available
Explosive limits : No data available

## 9.2. Other information

No additional information available

## SECTION 10: Stability and reactivity

#### 10.1. Reactivity

No additional information available

## 10.2. Chemical stability

Stable under normal conditions of storage, handling and use.

## 10.3. Possibility of hazardous reactions

This product is not hazardous.

### 10.4. Conditions to avoid

No additional information available

## 10.5. Incompatible materials

Strong oxidizing agent.

21/09/2017 Duchefa Biochemie B.V.

Safety Data Sheet

according to Regulation (EU) 2015/830

## 10.6. Hazardous decomposition products

When heated to decomposition, emits dangerous fumes. COx. NOx.

## **SECTION 11: Toxicological information**

## 11.1. Information on toxicological effects

Acute toxicity : Not classified

Skin corrosion/irritation : Not classified

Serious eye damage/irritation : Not classified

Respiratory or skin sensitisation : Not classified

Germ cell mutagenicity : Not classified

Carcinogenicity : Not classified

Reproductive toxicity : Not classified STOT-single exposure : Not classified STOT-repeated exposure : Not classified

Aspiration hazard : Not classified

## **SECTION 12: Ecological information**

### 12.1. Toxicity

No additional information available

## 12.2. Persistence and degradability

No additional information available

## 12.3. Bioaccumulative potential

DL-Dihydrozeatin (14894-18-9)	
Log Pow	0,466

## 12.4. Mobility in soil

No additional information available

## 12.5. Results of PBT and vPvB assessment

No additional information available

## 12.6. Other adverse effects

Additional information : Prevent entry to sewers and public waters. Avoid release to the environment

## **SECTION 13: Disposal considerations**

#### 13.1. Waste treatment methods

Waste treatment methods : Dispose in a safe manner in accordance with local/national regulations.

Avoid release to the environment.

21/09/2017 Duchefa Biochemie B.V. EN (English) SDS Ref.: D0933 5/7

Safety Data Sheet

according to Regulation (EU) 2015/830

## **SECTION 14: Transport information**

In accordance with ADR / IATA / IMDG

ADR	IMDG	IATA	
14.1. UN number			
Not applicable	Not applicable	Not applicable	
14.2. UN proper shipping name	14.2. UN proper shipping name		
Not applicable	Not applicable	Not applicable	
Not applicable	Not applicable	Not applicable	
14.3. Transport hazard class(es)			
Not applicable	Not applicable	Not applicable	
Not applicable	Not applicable	Not applicable	
14.4. Packing group			
Not applicable	Not applicable	Not applicable	
14.5. Environmental hazards			
Not applicable	Not applicable	Not applicable	
No supplementary information available			

#### 14.6. Special precautions for user

## . Overland transport

Not applicable

#### - Transport by sea

Not applicable

## - Air transport

Not applicable

## 14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

## **SECTION 15: Regulatory information**

## 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

## 15.1.1. EU-Regulations

No REACH Annex XVII restrictions

DL-Dihydrozeatin is not on the REACH Candidate List

DL-Dihydrozeatin is not on the REACH Annex XIV List

#### 15.1.2. National regulations

## Germany

12th Ordinance Implementing the Federal Immission Control Act -

: Is not subject of the 12. BlmSchV (Hazardous Incident Ordinance)

12.BImSchV

## **Netherlands**

SZW-lijst van kankerverwekkende

stoffen

: The substance is not listed

SZW-lijst van mutagene stoffen : The substance is not listed

 21/09/2017 Duchefa Biochemie B.V.
 EN (English)
 SDS Ref.: D0933
 6/7

## Safety Data Sheet

according to Regulation (EU) 2015/830

NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Borstvoeding : The substance is not listed

NIET-limitatieve lijst van voor de voortplanting giftige stoffen –

Voortplanting giftige s

: The substance is not listed

NIET-limitatieve lijst van voor de

voortplanting giftige stoffen – Ontwikkeling

: The substance is not listed

## 15.2. Chemical safety assessment

No additional information available

## **SECTION 16: Other information**

## Indication of changes:

1.4	Emergency number	Modified	
8.2	Hand protection	Modified	Specified material, thickness, et
			cetera of gloves

#### Abbreviations and acronyms:

	<b>,</b> -
ATE	Acute Toxicity Estimate
ADR	European Agreement concerning the International Carriage of Dangerous Goods by Road
BCF	Bioconcentration factor
CLP	Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008
DPD	Dangerous Preparations Directive 1999/45/EC
DSD	Dangerous Substances Directive 67/548/EEC
IATA	International Air Transport Association
IMDG	International Maritime Dangerous Goods
LC50	Median lethal concentration
LD50	Median lethal dose
LOAEL	Lowest Observed Adverse Effect Level
NOAEC	No-Observed Adverse Effect Concentration
PBT	Persistent Bioaccumulative Toxic
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006
SDS	Safety Data Sheet

#### Data sources

: Manufacturer. REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006.

## SDS Biochemicals version 2017

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product

21/09/2017 Duchefa Biochemie B.V. EN (English) SDS Ref.: D0933 7/7