

Safety Data Sheet

according to Regulation (EU) 2015/830

Version: 2.0

according to Regulation (EU) 2015/830

Date of issue: 23/03/2012 Revision date 15/12/2017

Supersedes 04/04/2012

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### 1.1. Product identifier

Product form : Article

Trade name : 6-y-y-(Dimethylallylamino)purine riboside

EC-No. : 231-771-3 CAS-No. : 7724-76-7 Product code : D0934

Formula : C15H21N5O4

Synonyms : N6-(2-Isopentenyl)adenosine / 2-iP-riboside

#### 1.2. Relevant identified uses of the substance or mixture and uses advised against

#### 1.2.1. Relevant identified uses

Main use category : Professional use

Industrial/Professional use spec : For professional use only. Duchefa Biochemie B.V. products are intended

only for "in vitro laboratory" research purposes.

#### 1.2.2. Uses advised against

No additional information available

#### 1.3. Details of the supplier of the safety data sheet

#### Supplier

Duchefa Biochemie B.V.
A. Hofmanweg 71
2031 BH Haarlem - The Netherlands
T +31(0)23-5319093 - F +31(0)23-5318027
info@duchefa.nl

#### 1.4. Emergency telephone number

Emergency number : Supplier contact information:

+31(0)23-5319093 (M-F 09:00-17:00) +31(0)6-30109355 (outside office hours)

SDS Ref.: D0934

Organisation/Company	Address	Comment
World Health Organization world	http://apps.who.int/poisoncentres/	Consult website for a local poison centre
directory of poison centres		·

#### **SECTION 2: Hazards identification**

#### 2.1. Classification of the substance or mixture

### Classification according to Regulation (EC) No. 1272/2008 [CLP]

Not classified

#### Adverse physicochemical, human health and environmental effects

No additional information available

Safety Data Sheet

according to Regulation (EU) 2015/830

#### 2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

No labelling applicable

#### 2.3. Other hazards

No additional information available

### **SECTION 3: Composition/information on ingredients**

#### 3.1. Substances

Not applicable

#### 3.2. Mixtures

Name	Product identifier	%	Classification according to Regulation (EC) No. 1272/2008 [CLP]
6-y-y-(Dimethylallylamino)purine riboside	(CAS-No.) 7724-76-7 (EC-No.) 231-771-3	> 97	Not classified

### **SECTION 4: First aid measures**

#### 4.1. Description of first aid measures

First-aid measures general : Seek medical attention if ill effect or irritation develops.

First-aid measures after inhalation : Move to fresh air.

First-aid measures after skin contact : Wash skin with mild soap and water.

First-aid measures after eye contact : Rinse with water. First-aid measures after ingestion : Rinse mouth.

#### 4.2. Most important symptoms and effects, both acute and delayed

No additional information available

#### 4.3. Indication of any immediate medical attention and special treatment needed

No additional information available

# SECTION 5: Firefighting measures

#### 5.1. Extinguishing media

Suitable extinguishing media : Alcohol resistant foam

Dry chemical powder Carbon dioxide (CO2)

Water spray.

### 5.2. Special hazards arising from the substance or mixture

Hazardous decomposition products in : Under fire conditions, hazardous fumes will be present. COx. NOx.

case of fire

#### 5.3. Advice for firefighters

Protection during firefighting : Wear proper protective equipment.

Other information : Prevent fire fighting water from entering the environment.

 15/12/2017 Duchefa Biochemie B.V.
 EN (English)
 SDS Ref.: D0934
 2/7

Safety Data Sheet

according to Regulation (EU) 2015/830

### SECTION 6: Accidental release measures

#### 6.1. Personal precautions, protective equipment and emergency procedures

#### 6.1.1. For non-emergency personnel

Protective equipment : Wear suitable protective clothing.

Measures in case of dust release : Avoid dust formation.

#### 6.1.2. For emergency responders

No additional information available

#### 6.2. Environmental precautions

Prevent entry to sewers and public waters.

### 6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Sweep up dry powder and dispose properly.

#### 6.4. Reference to other sections

No additional information available

### **SECTION 7: Handling and storage**

#### 7.1. Precautions for safe handling

Precautions for safe handling : Store in dry, well-ventilated area.

#### 7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Keep container tightly closed and protected from light

Store in dry, cool, well-ventilated area

Keep cool.

Storage temperature : 2 - 8 °C

#### 7.3. Specific end use(s)

For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes.

### SECTION 8: Exposure controls/personal protection

#### 8.1. Control parameters

No additional information available

#### 8.2. Exposure controls

#### Hand protection:

Туре	Material	Permeation	Thickness (mm)	Standard
Gloves	Nitrile rubber (NBR)	6 (> 480 minutes)	0,11	EN 374

#### Eye protection:

Safety glasses (to European standard EN 166 or equivalent)

#### Skin and body protection:

Wear suitable protective clothing

#### Respiratory protection:

Where excessive dust may result, wear approved mask. Type P1 (EN 143)

Safety Data Sheet

according to Regulation (EU) 2015/830

### SECTION 9: Physical and chemical properties

#### 9.1. Information on basic physical and chemical properties

Physical state : Solid

Appearance : Powder. Crystals.

Molecular mass : 335,4 g/mol

Colour : White.

Odour : Negligible.

Odour threshold : No data available pH : No data available Relative evaporation rate : No data available

(butylacetate=1)

Melting point : 140 - 150 °C Freezing point : No data available : No data available **Boiling point** : No data available Flash point Auto-ignition temperature : No data available Decomposition temperature : No data available : No data available Flammability (solid, gas) : No data available Vapour pressure Relative vapour density at 20 °C : No data available Relative density : No data available

Solubility : Readily soluble in water.

Log Pow : No data available
Viscosity, kinematic : No data available
Viscosity, dynamic : No data available
Explosive properties : No data available
Oxidising properties : No data available
Explosive limits : No data available

#### 9.2. Other information

No additional information available

### SECTION 10: Stability and reactivity

#### 10.1. Reactivity

No additional information available

#### 10.2. Chemical stability

No additional information available

#### 10.3. Possibility of hazardous reactions

No additional information available

#### 10.4. Conditions to avoid

Light (daylight).

#### 10.5. Incompatible materials

Strong oxidizers.

#### 10.6. Hazardous decomposition products

When heated to decomposition, emits dangerous fumes.

COx

Safety Data Sheet

according to Regulation (EU) 2015/830

NOx.

## **SECTION 11: Toxicological information**

#### 11.1. Information on toxicological effects

Acute toxicity : Not classified

Skin corrosion/irritation : Not classified

Serious eye damage/irritation : Not classified

Respiratory or skin sensitisation : Not classified

Germ cell mutagenicity : Not classified

Carcinogenicity : Not classified

Reproductive toxicity : Not classified STOT-single exposure : Not classified STOT-repeated exposure : Not classified

Aspiration hazard : Not classified

### SECTION 12: Ecological information

#### 12.1. Toxicity

No additional information available

#### 12.2. Persistence and degradability

No additional information available

#### 12.3. Bioaccumulative potential

No additional information available

#### 12.4. Mobility in soil

No additional information available

#### 12.5. Results of PBT and vPvB assessment

No additional information available

#### 12.6. Other adverse effects

No additional information available

### **SECTION 13: Disposal considerations**

#### 13.1. Waste treatment methods

Waste treatment methods : Dispose in a safe manner in accordance with local/national regulations.

### **SECTION 14: Transport information**

In accordance with ADN / ADR / IATA

ADR	IATA	ADN
14.1. UN number		
Not regulated	Not regulated	Not applicable
14.2. UN proper shipping name		
Not regulated	Not regulated	Not applicable

Safety Data Sheet

according to Regulation (EU) 2015/830

ADR	IATA	ADN
Not regulated	Not regulated	Not applicable
14.3. Transport hazard class(es)		
Not regulated	Not regulated	Not applicable
Not regulated	Not regulated	Not applicable
14.4. Packing group		
Not regulated	Not regulated	Not applicable
14.5. Environmental hazards		
Not regulated	Not regulated	Not applicable
No supplementary information available		

#### 14.6. Special precautions for user

#### Overland transport

Not regulated

#### - Air transport

Not regulated

#### - Inland waterway transport

Not applicable

#### Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

### **SECTION 15: Regulatory information**

#### Safety, health and environmental regulations/legislation specific for the substance or mixture

#### 15.1.1. EU-Regulations

Contains no substance on the REACH candidate list

#### 15.1.2. National regulations

Ensure all national/local regulations are observed.

#### Germany

12th Ordinance Implementing the Federal Immission Control Act -: Is not subject of the 12. BlmSchV 12.BImSchV (Hazardous Incident Ordinance)

#### **Netherlands**

SZW-lijst van kankerverwekkende stoffen : None of the components are listed SZW-lijst van mutagene stoffen : None of the components are listed NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Borstvoeding : None of the components are listed NIET-limitatieve lijst van voor de voortplanting giftige stoffen -: None of the components are listed Vruchtbaarheid

#### NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Ontwikkeling : None of the components are listed

#### Chemical safety assessment

No additional information available

#### **SECTION 16: Other information**

Indication of changes:

Safety Data Sheet

according to Regulation (EU) 2015/830

1.4	Emergency number	Modified	
8.2	Hand protection	Modified	Specified material, thickness, et cetera of gloves
9.1	Melting point	Modified	

#### Abbreviations and acronyms:

ATE	Acute Toxicity Estimate
ADR	European Agreement concerning the International Carriage of Dangerous Goods by Road
BCF	Bioconcentration factor
CLP	Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008
DPD	Dangerous Preparations Directive 1999/45/EC
DSD	Dangerous Substances Directive 67/548/EEC
IATA	International Air Transport Association
IMDG	International Maritime Dangerous Goods
LC50	Median lethal concentration
LD50	Median lethal dose
LOAEL	Lowest Observed Adverse Effect Level
NOAEC	No-Observed Adverse Effect Concentration
PBT	Persistent Bioaccumulative Toxic
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006
SDS	Safety Data Sheet

Data sources

: Manufacturer. ECHA (European Chemicals Agency).

#### SDS Biochemicals version 2017

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product