

## FICHE TECHNIQUE

### **D1342: DEXTRAN SULPHATE SODIUM**

CAS number 9011-18-1

Molecular weight  $(C_6H_9O_5SO_3Na)_{n+1}$

Free sulphate < 0.2 %

pH aqueous solution (1 %) 5.0 - 7.5

Clarity (15 % solution) : No suspended particles

#### **Storage condition**

- Store at room temperature
- Soluble in hot water

D1342.0010/ 10 g

D1342.0050/ 50 g

D1342.0100/ 100 g