

FICHE TECHNIQUE

E0122 : ERYTHROMYCIN



CAS number 114-07-8
Molecular weight C₃₇H₆₇NO₁₃ = 733.9
Assay > 93 %
Water < 6.5 %

Storage condition

- Store dry at room temperature
- Soluble in ethanol

Hazard statements

- **H317** May cause an allergic skin reaction
- **H334** May cause allergy or asthma symptoms or breathing difficulties if inhaled

Precautionary statements

- **P261** Avoid breathing dust/fume/gas/mist/vapours/spray
- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P311** Call a POISON CENTER or doctor/physician
- **P342** If experiencing respiratory symptoms:

E0122.0010/ 10 g
E0122.0025/ 25 g