

FICHE TECHNIQUE



E0207 : ERIKSSON (ER) MEDIUM

Micro Elements

mg/l

µM

CoCl₂.6H₂O

0.0025

0.01

CuSO₄.5H₂O

0.0025

0.01

FeNaEDTA

36.70

100.00

H₃BO₃

0.63

10.19

MnSO₄.H₂O

1.69

10.00

Na₂MoO₄.2H₂O

0.025

0.10

ZnSO₄.7H₂O

9.15

31.80

Macro Elements

mg/l

mM

CaCl₂

332.02

2.99

KH₂PO₄

340.00

2.50

KNO₃

1900.00

18.79

MgSO₄

180.54

1.50

NH₄NO₃

1200.00

14.99

Eriksson T., *Physiol. Plant*, 18, 976 (1965).

Hazard statements

- **H315** Causes skin irritation
- **H319** Causes serious eye irritation
- **H335** May cause respiratory irritation

Precautionary statements

- **P261** Avoid breathing dust/fume/gas/mist/vapours/spray
- **P302** IF ON SKIN:
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes
- **P352** Wash with soap and water

E0207.0001/ 1 l

E0207.0010/ 10 l

E0207.0050/ 50 l

Parc Euromédecine -20, rue Robert Koch - BP 27253 34086 MONTPELLIER Cedex 4 - Tel : +33(0)4 67 10 59 00

RCS 461 800 559 B Montpellier - SIRET 461 800 559 00035 - CODE T.V.A CEE. : FR 62 461 800 559

contact@groupe-meridis.fr - www.groupe-meridis.fr