

FICHE TECHNIQUE



E0208 : ERIKSSON (ER) MEDIUM INCLUDING VITAMINS

Micro Elements

mg/l

μM

CoCl₂.6H₂O

0.0025

0.01

CuSO₄.5H₂O

0.0025

0.01

FeNaEDTA

36.70

100.00

H₃BO₃

0.63

10.19

MnSO₄.H₂O

1.69

10.00

Na₂MoO₄.2H₂O

0.025

0.10

ZnSO₄.7H₂O

9.15

31.80

Macro Elements

mg/l

mM

CaCl₂

332.02

2.99

KH₂PO₄

340.00
2.50
KNO3
1900.00
18.79
MgSO4
180.54
1.50
NH4NO3
1200.00
14.99

Vitamins

mg/l

µM

Glycine

2.00

26.64

Nicotinic acid

0.50

4.06

Pyridoxine HCl

0.50

2.43

Thiamine HCl

0.50

1.48

Total concentration Micro and Macro elements including -vitamins: 4004.42 mg/l

Eriksson T., Physiol. Plant, 18, 976 (1965).

Hazard statements

- **H315** Causes skin irritation
- **H319** Causes serious eye irritation
- **H335** May cause respiratory irritation

Precautionary statements

- **P261** Avoid breathing dust/fume/gas/mist/vapours/spray
- **P302** IF ON SKIN:
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes
- **P352** Wash with soap and water

E0208.0001/ 1 l

E0208.0010/ 10 l

E0208.0050/ 50 l