

FICHE TECHNIQUE



E0402: ERIKSSON (ER) VITAMIN MIXTURE

Vitamins

mg/l

µM

Glycine

2.00

26.64

Nicotinic acid

0.50

4.06

Pyridoxine HCl

0.50

2.43

Thiamine HCl

0.50

1.48

Eriksson T., *Physiol. Plant*, 18, 976 (1965).

Hazard statements

- **H319** Causes serious eye irritation

Precautionary statements

- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes

E0402.0100/ >100 ml